
DICOMweb Client Documentation

Release 0.59.1

Markus D. Herrmann

Feb 27, 2023

CONTENTS:

1	Introduction	1
2	Installation guide	3
2.1	Requirements	3
2.2	Installation	3
3	User guide	5
3.1	Application Programming Interface (API)	5
3.2	Command Line Interface (CLI)	14
4	Developer guide	17
4.1	Pull requests	17
4.2	Coding style	17
4.3	Running tests	17
4.4	Building documentation	18
5	Conformance statement	19
5.1	QIDO-RS	19
5.2	WADO-RS	19
5.3	STOW-RS	19
6	License	21
7	API documentation	23
7.1	dicomweb_client package	23
8	Indices and tables	71
	Python Module Index	73
	Index	75

INTRODUCTION

The `dicomweb-client` build distribution provides client interfaces for DICOMweb RESTful services QIDO-RS, WADO-RS and STOW-RS to search, retrieve and store DICOM objects over the web, respectively. For more information about DICOMweb please refer to the documentation of the [DICOM standard](#), in particular [PS3.18](#).

The `dicomweb_client` Python package exposes

- Application Programming Interface (API) (see [api](#) module)
- Command Line Interface (CLI) (see [cli](#) module)

INSTALLATION GUIDE

2.1 Requirements

- Python (version 3.5 or higher)
- Python package manager `pip`

For support of image formats:

- JPEG (`libjpeg` or `libjpeg-turbo`)
- JPEG2000 (`openjpeg`)
- PNG (`libpng`)

2.2 Installation

Pre-build package available at PyPi:

```
pip install dicomweb-client
```

Additional dependencies required for extensions compatible with [Google Cloud Platform \(GCP\)](#) may be installed as:

```
pip install dicomweb-client[gcp]
```

Source code available at Github:

```
git clone https://github.com/ImagingDataCommons/dicomweb-client ~/dicomweb-client  
pip install ~/dicomweb-client
```


The client can be used with any DICOMweb server, such as [dcm4che](#), [orthanc](#) or [DICOMcloud](#).

3.1 Application Programming Interface (API)

3.1.1 Interacting with a remote DICOMweb server

To interact with a publicly accessible server, you only need to provide the url for the server address.

```
from dicomweb_client.api import DICOMwebClient

client = DICOMwebClient(url="https://mydicomwebserver.com")
```

Some servers expose the different DICOMweb RESTful services using different path prefixes. For example, the publicly accessible [DICOMcloud server](#) uses the prefixes "qidors", "wadors", and "stowrs" for QIDO-RS, WADO-RS, and STOW-RS, respectively. You can specify these prefixes using `qido_url_prefix`, `wado_url_prefix`, and `stow_url_prefix`.

```
from dicomweb_client.api import DICOMwebClient

client = DICOMwebClient(
    url="https://dicomcloud.azurewebsites.net",
    qido_url_prefix="qidors",
    wado_url_prefix="wadors",
    stow_url_prefix="stowrs"
)
```

3.1.2 Authentication and authorization

To interact with servers requiring authentication, `DICOMwebClient` accepts arbitrary authentication handlers derived from `requests.auth.AuthBase` (see [here](#) for details).

```
from requests.auth import HTTPBasicAuth
from dicomweb_client.api import DICOMwebClient
from dicomweb_client.session_utils import create_session_from_auth

auth = HTTPBasicAuth('myusername', 'mypassword')
session = create_session_from_auth(auth)
```

(continues on next page)

(continued from previous page)

```
client = DICOMwebClient(
    url="https://mydicomwebserver.com",
    session=session
)
```

To simplify usage for basic HTTP authentication, you may also directly provide a username and password using the corresponding arguments.

```
from dicomweb_client.api import DICOMwebClient
from dicomweb_client.session_utils import create_session_from_user_pass

session = create_session_from_user_pass(
    username='myusername',
    password='mypassword'
)

client = DICOMwebClient(
    url="https://mydicomwebserver.com",
    session=session
)
```

To interact with servers supporting token-based authorization, you can provide the access token using the headers argument (the header will be included in every client request message).

```
from dicomweb_client.api import DICOMwebClient

access_token = "mytoken"
client = DICOMwebClient(
    url="https://mydicomwebserver.com",
    headers={"Authorization": "Bearer {}".format(access_token)}
)
```

To interact with servers requiring certificate-based authentication, you can provide the CA bundle and client certificate using the ca_bundle and cert arguments, respectively.

```
from dicomweb_client.api import DICOMwebClient
from dicomweb_client.session_utils import (
    create_session,
    add_certs_to_session
)

session = create_session()
session = add_certs_to_session(
    session=session,
    ca_bundle="/path/to/ca.crt",
    cert="/path/to/cert.pem"
)

client = DICOMwebClient(url="https://mydicomwebserver.com")
```

To interact with a server of the Google Healthcare API requiring OpenID Connect based authentication and authorization, provide a session authenticated using the Google Cloud Platform (GCP) credentials. See [GCP documentation](#) for

details.

The library provides the `gcp` extension, which facilitates interacting with the DICOMweb interface of the Google Healthcare API. Note that the `gcp` extension is optional and requires installation of the package distribution with the `gcp` extra requirements: `$ pip install dicomweb-client[gcp]`.

```
from dicomweb_client.api import DICOMwebClient
from dicomweb_client.ext.gcp.session_utils import create_session_from_gcp_credentials
from dicomweb_client.ext.gcp.uri import GoogleCloudHealthcareURL

session = create_session_from_gcp_credentials()

url = GoogleCloudHealthcareURL(
    project_id='my-project',
    location='us-east4',
    dataset_id='my-dataset',
    dicom_store_id='my-store'
)

client = DICOMwebClient(
    url=str(url),
    session=session
)
```

3.1.3 Accessing local DICOM Part10 files

The package provides the `dicomweb_client.api.DICOMfileClient` class for accessing data stored as DICOM Part10 files on a file system. The class exposes the same `dicomweb_client.api.DICOMClient` interface as the `dicomweb_client.api.DICOMwebClient` and can be used as a drop-in replacement.

To create a `dicomweb_client.api.DICOMfileClient` instance, you should pass the location of the desired data store on the file system. Note that the `file://` scheme designator is *required*.

```
from dicomweb_client.api import DICOMfileClient

client = DICOMfileClient("file:///path/to/directory")
```

3.1.4 STOW-RS StoreInstances

Store a single dataset obtained from a PS3.10 file:

```
import pydicom

filename = "/path/to/file.dcm"
dataset = pydicom.dcmread(filename)

client.store_instances(datasets=[dataset])
```

3.1.5 QIDO-RS SearchForStudies

Search for all studies (up to server-defined maximum set per call - see below to iteratively get all studies):

```
studies = client.search_for_studies()
```

Search for studies filtering by *PatientID*:

```
studies = client.search_for_studies(search_filters={'PatientID': 'ABC123'})
```

Note that attributes can be specified in `search_filters` using either the keyword or the tag:

```
studies = client.search_for_studies(search_filters={'00100020': 'ABC123'})
```

Search for all studies but limit the number of returned results using the `limit` parameter.

```
studies_subset = client.search_for_studies(limit=100)
```

A server may also automatically limit the number of results that it returns per search request. In this case, the method can be called repeatedly to request remaining results using the `offset` parameter.

```
studies = []
offset = 0
while True:
    subset = client.search_for_studies(offset=offset)
    if len(subset) == 0:
        break
    studies.extend(subset)
    offset += len(subset)
```

The same can be achieved more conveniently using the `get_remaining` parameter.

```
studies = client.search_for_studies(get_remaining=True)
```

3.1.6 QIDO-RS SearchForSeries

Search for all series:

```
series = client.search_for_series()
```

Search for series of a given study:

```
series = client.search_for_series('1.2.826.0.1.3680043.8.1055.1.20111103111148288.
↪98361414.79379639')
```

Search for series filtering by *AccessionNumber*:

```
series = client.search_for_series(search_filters={'AccessionNumber': '123456'})
```

Search for series filtering by *AccessionNumber* (using wildcard `?` to match a range of numbers):

```
series = client.search_for_series(search_filters={'AccessionNumber': '12345?'})
```

Search for series filtering by *SeriesDescription*:

```
series = client.search_for_series(search_filters={'SeriesDescription': 'T2 AXIAL'})
```

Search for series filtering by *SeriesDescription* (using wildcard * to match a range of descriptions):

```
series = client.search_for_series(search_filters={'SeriesDescription': 'T2 AX*'})
```

Search for series filtering by *Modality*:

```
series = client.search_for_series(search_filters={'Modality': 'SM'})
```

3.1.7 QIDO-RS SearchForInstances

Search for all instances:

```
instances = client.search_for_instances()
```

Search for instances of a given study and series:

```
instances = client.search_for_instances(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↪',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↪'
)
```

Search for instances filtering by *SOPClassUID*:

```
instances = client.search_for_instances(search_filters={'SOPClassUID': '1.2.840.10008.5.
↪1.4.1.1.2'})
```

3.1.8 WADO-RS RetrieveStudy

Retrieve instances of a given study:

```
instances = client.retrieve_study('1.2.826.0.1.3680043.8.1055.1.20111103111148288.
↪98361414.79379639')
```

3.1.9 WADO-RS RetrieveSeries

Retrieve instances of a given series:

```
instances = client.retrieve_series(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↪',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↪'
)
```

Retrieve full instances of a given series using specific JPEG 2000 transfer syntax for encoding of bulk data:

```
instance = client.retrieve_instance(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↔',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↔',
    media_types=(('application/dicom', '1.2.840.10008.1.2.4.90', ), )
)
```

Retrieve bulk data of instances of a given series using specific JPEG 2000 transfer syntax:

```
instance = client.retrieve_instance(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↔',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↔',
    media_types=(('image/jp2', '1.2.840.10008.1.2.4.90', ), )
)
```

3.1.10 WADO-RS RetrieveInstance

Retrieve full instance using default Explicit VR Little Endian transfer syntax for encoding of bulk data:

```
instance = client.retrieve_instance(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↔',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↔',
    sop_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534'
)
```

Retrieve full instance using specific JPEG 2000 transfer syntax for encoding of bulk data:

```
instance = client.retrieve_instance(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↔',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↔',
    sop_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534',
    media_types=(('application/dicom', '1.2.840.10008.1.2.4.90', ), )
)
```

Retrieve bulk data of instance using specific JPEG 2000 transfer syntax:

```
instance = client.retrieve_instance(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↔',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↔',
    sop_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534',
    media_types=(('image/jp2', '1.2.840.10008.1.2.4.90', ), )
)
```

3.1.11 WADO-RS RetrieveMetadata

Retrieve metadata for instances of a given study:

```
metadata = client.retrieve_study_metadata('1.2.826.0.1.3680043.8.1055.1.
↳20111103111148288.98361414.79379639')
```

Retrieve metadata for instances of a given series:

```
metadata = client.retrieve_series_metadata(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↳',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↳'
)
```

Retrieve metadata for a particular instance:

```
metadata = client.retrieve_instance_metadata(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↳',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↳',
    sop_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534'
)
```

Note: WADO-RS RetrieveMetadata always returns metadata at the instance-level, `retrieve_study_metadata()` and `retrieve_series_metadata()` return an array of metadata items for each instance belonging to a given study and series, respectively.

3.1.12 WADO-RS RetrieveFrames

Retrieve a set of frames with default transfer syntax (“application/octet-stream”):

```
frames = client.retrieve_instance_frames(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↳',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↳',
    sop_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534',
    frame_numbers=[1, 2]
)
```

Retrieve a set of frames of a given instances as JPEG compressed image:

```
frames = client.retrieve_instance_frames(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↳',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↳',
    sop_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534',
```

(continues on next page)

(continued from previous page)

```

frame_numbers=[1, 2],
media_types=('image/jpeg', )
)

```

Retrieve a set of frames of a given instances as compressed image in any available format:

```

frames = client.retrieve_instance_frames(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↔',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↔',
    sop_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534',
    frame_numbers=[1, 2],
    media_types=('image/*', )
)

```

Retrieve a set of frames of a given instances as either JPEG 2000 or JPEG-LS compressed image:

```

frames = client.retrieve_instance_frames(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↔',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↔',
    sop_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534',
    frame_numbers=[1, 2],
    media_types=('image/jp2', 'image/jls', )
)

```

Retrieve a set of frames of a given instances as either JPEG, JPEG 2000 or JPEG-LS lossless compressed image using specific transfer syntaxes:

```

frames = client.retrieve_instance_frames(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
↔',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
↔',
    sop_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534',
    frame_numbers=[1, 2],
    media_types=(
        ('image/jpeg', '1.2.840.10008.1.2.4.57', ),
        ('image/jp2', '1.2.840.10008.1.2.4.90', ),
        ('image/jls', '1.2.840.10008.1.2.4.80', ),
    )
)

```


3.1.13 WADO-RS RetrieveBulkdata

Retrieve bulk data given a URL:

```
data = client.retrieve_bulkdata('https://mydicomwebserver.com/studies/...')
```

3.1.14 WADO-RS RetrieveRenderedTransaction

Retrieve a single-frame image instance rendered as a PNG compressed image:

```
frames = client.retrieve_instance_rendered(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
    ↪',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
    ↪',
    sop_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534',
    media_types=('image/png', )
)
```

Retrieve a single frame of a multi-frame image instance rendered as a high-quality JPEG compressed image that includes an ICC profile:

```
frames = client.retrieve_instance_frames_rendered(
    study_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639
    ↪',
    series_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034
    ↪',
    sop_instance_uid='1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534',
    frame_numbers=[1],
    media_types=('image/jpeg', ),
    params={'quality': 95, 'iccprofile': 'yes'}
)
```

When frames are retrieved in image format, they can be converted into a *NumPy* array using the *PIL* module:

```
from io import BytesIO

import numpy as np
from PIL import Image

image = Image.open(BytesIO(frames[0]))
array = np.array(image)
```

Warning: Retrieving images using lossy compression methods may lead to image recompression artifacts if the images have been stored lossy compressed.

3.1.15 Loading JSON Data To pydicom

Load metadata from JSON format into a `pydicom.dataset.Dataset` object. A common use for this is translating metadata received from a `RetrieveMetadata` or a `SearchFor`-style request:

```
from pydicom.dataset import Dataset

metadata = client.retrieve_study_metadata('1.2.826.0.1.3680043.8.1055.1.
↪20111103111148288.98361414.79379639')
metadata_datasets = [Dataset.from_json(ds) for ds in metadata]
```

Note that the metadata may include references to `BulkData` elements. By default, `BulkData` elements will not be handled and the values not be automatically retrieved. To handle `BulkData` elements and retrieve their values, one has to provide a `bulk_data_uri_handler` callable to the `pydicom.dataset.Dataset.from_json()` method.

3.2 Command Line Interface (CLI)

Search for studies:

```
dicomweb_client --url https://dicomcloud.azurewebsites.net/qidors search studies
```

Retrieve metadata for all instances of a given study:

```
dicomweb_client --url https://dicomcloud.azurewebsites.net/wadors \
  retrieve studies \
  --study 1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639 \
  metadata
```

The output can be *dicomized* for human interpretation:

```
dicomweb_client --url https://dicomcloud.azurewebsites.net/wadors \
  retrieve studies \
  --study 1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639 \
  metadata \
  --dicomize
```

Retrieve the full Part 3.10 files for all instances of a given study:

```
dicomweb_client --url https://dicomcloud.azurewebsites.net/wadors \
  retrieve studies \
  --study 1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639 \
  full
```

Retrieve a single frame of a given instances as JPEG compressed image:

```
dicomweb_client --url https://dicomcloud.azurewebsites.net/wadors \
  retrieve instances \
  --study 1.2.826.0.1.3680043.8.1055.1.20111103111148288.98361414.79379639 \
  --series 1.2.826.0.1.3680043.8.1055.1.20111103111208937.49685336.24517034 \
  --instance 1.2.826.0.1.3680043.8.1055.1.20111103111208937.40440871.13152534 \
  frames \
  --numbers 1 \
  --media-type image/jpeg
```

Store instances to a Google DICOMweb store:

```
dicomweb_client --url https://healthcare.googleapis.com/v1beta1/projects/MYPROJECT/  
↪locations/us-central1/datasets/MYDATASET/dicomStores/MYDICOMSTORE/dicomWeb \  
  --token $(gcloud auth print-access-token) \  
  store instances \  
  dicomfiles/*
```


DEVELOPER GUIDE

Source code is available at Github and can be cloned via git:

```
git clone https://github.com/ImagingDataCommons/dicomweb-client ~/dicomweb-client
```

The `dicomweb_client` package can be installed in *develop* mode for local development:

```
pip install -e ~/dicomweb-client
```

4.1 Pull requests

Don't commit code changes to the `master` branch. New features should be implemented in a separate branch called `feature/*` and bug fixes should be applied in separate branch called `bugfix/*`.

Before creating a pull request on Github, read the coding style guideline, run the tests and check PEP8 compliance.

4.2 Coding style

Code must comply with [PEP 8](#). The `flake8` package is used to enforce compliance.

The project uses `numpydoc` for documenting code according to [PEP 257](#) docstring conventions. Further information and examples for the NumPy style can be found at the [NumPy Github repository](#) and the website of the [Napoleon sphinx extension](#).

All API classes, functions and modules must be documented (including "private" functions and methods). Each docstring must describe input parameters and return values. Types must be specified using type hints as specified by [PEP 484](#) (see `typing` module) in both the function definition as well as the docstring.

4.3 Running tests

The project uses `pytest` to write and runs unit tests. Tests should be placed in a separate `tests` folder within the package root folder. Files containing actual test code should follow the pattern `test_*.py`.

Install requirements:

```
pip install -r ~/dicomweb-client/requirements_test.txt
```

Run tests (including checks for PEP8 compliance):

```
cd ~/dicomweb-client
pytest --flake8
```

4.4 Building documentation

Install requirements:

```
pip install -r ~/dicomweb-client/requirements_docs.txt
```

Build documentation in *HTML* format:

```
cd ~/dicomweb-client
sphinx-build -b html docs/ docs/build/
```

The built `index.html` file will be located in `docs/build`.

CONFORMANCE STATEMENT**5.1 QIDO-RS**

Method	Resource	Implemented
GET	SearchForStudies	Y
GET	SearchForSeries	Y
GET	SearchForInstances	Y

5.2 WADO-RS

Method	Resource	Implemented
GET	RetrieveStudy	Y
GET	RetrieveSeries	Y
GET	RetrieveInstance	Y
GET	RetrieveMetadata	Y*
GET	RetrieveBulkdata	Y
GET	RetrieveFrames	Y
GET	RetrieveRenderedTransaction	Y

- *Metadata* resource representations are requested in JSON format according to the [DICOM JSON model](#) using `application/dicom+json` media type. Retrieval of metadata in XML form using `application/dicom+xml` is not supported.

5.3 STOW-RS

Method	Resource	Implemented
POST	StoreInstances	Y

LICENSE

DICOMweb Client is free and open source software licensed under the permissive [MIT license](#).

7.1 dicomweb_client package

7.1.1 dicomweb_client.api module

Application Programming Interface (API)

class `dicomweb_client.api.DICOMClient(*args, **kwargs)`

Bases: `Protocol`

Protocol for DICOM clients based on DICOMweb interface.

base_url: `str`

delete_instance(*study_instance_uid*, *series_instance_uid*, *sop_instance_uid*)

Delete specified instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID

Returns

HTTP response object returned.

Return type

`requests.models.Response`

Note: The Delete Instance resource is not part of the DICOM standard and may not be supported by all origin servers.

delete_series(*study_instance_uid*, *series_instance_uid*)

Delete all instances of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID

Returns

HTTP response object returned.

Return typerequests.models.Response

Note: The Delete Series resource is not part of the DICOM standard and may not be supported by all origin servers.

delete_study(*study_instance_uid*)

Delete all instances of a study.

Parameters**study_instance_uid** (*str*) – Study Instance UID**Returns**

HTTP response object returned.

Return typerequests.models.Response

Note: The Delete Study resource is not part of the DICOM standard and may not be supported by all origin servers.

delete_url_prefix: Optional[*str*] = None**iter_bulkdata**(*url*, *media_types=None*, *byte_range=None*)

Iterate over bulk data items at a given location.

Parameters

- **url** (*str*) – Location of the bulk data
- **media_types** (*Union[Tuple[Union[*str*, Tuple[*str*, *str*]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes
- **byte_range** (*Union[Tuple[int, int], None]*, *optional*) – Start and end of byte range

Returns

Bulk data items

Return type

Iterator[bytes]

iter_instance_frames(*study_instance_uid*, *series_instance_uid*, *sop_instance_uid*, *frame_numbers*, *media_types=None*)

Iterate over frames of an image instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **frame_numbers** (*Sequence[int]*) – One-based positional indices of the frames within the instance

- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, optional) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes

Returns

Pixel data for each frame

Return type

Iterator[bytes]

iter_series(*study_instance_uid, series_instance_uid, media_types=None*)

Iterate over all instances of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, optional) – Acceptable media types and optionally the UIDs of the acceptable transfer syntaxes

Returns

Instances

Return type

Iterator[pydicom.dataset.Dataset]

Note: Instances are by default retrieved using Implicit VR Little Endian transfer syntax (Transfer Syntax UID "1.2.840.10008.1.2"). This means that Pixel Data of Image instances will be retrieved uncompressed. To retrieve instances in any available transfer syntax (typically the one in which instances were originally stored), specify acceptable transfer syntaxes using the wildcard ("application/dicom", "*").

iter_study(*study_instance_uid, media_types=None*)

Iterate over all instances of a study.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, optional) – Acceptable media types and optionally the UIDs of the acceptable transfer syntaxes

Returns

Instances

Return type

Iterator[pydicom.dataset.Dataset]

Note: Instances are by default retrieved using Implicit VR Little Endian transfer syntax (Transfer Syntax UID "1.2.840.10008.1.2"). This means that Pixel Data of Image instances will be retrieved uncompressed. To retrieve instances in any available transfer syntax (typically the one in which instances were originally stored), specify acceptable transfer syntaxes using the wildcard ("application/dicom", "*").

protocol: Optional[str]

qido_url_prefix: Optional[str] = None

retrieve_bulkdata(url, media_types=None, byte_range=None)

Retrieve bulk data at a given location.

Parameters

- **url** (str) – Location of the bulk data
- **media_types** (Union[Tuple[Union[str, Tuple[str, str]], ...], None], optional) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes
- **byte_range** (Union[Tuple[int, int], None], optional) – Start and end of byte range

Returns

Bulk data items

Return type

List[bytes]

retrieve_instance(study_instance_uid, series_instance_uid, sop_instance_uid, media_types=None)

Retrieve an individual instance.

Parameters

- **study_instance_uid** (str) – Study Instance UID
- **series_instance_uid** (str) – Series Instance UID
- **sop_instance_uid** (str) – SOP Instance UID
- **media_types** (Union[Tuple[Union[str, Tuple[str, str]], ...], None], optional) – Acceptable media types and optionally the UIDs of the acceptable transfer syntaxes

Returns

Instance

Return type

pydicom.dataset.Dataset

Note: Instances are by default retrieved using Implicit VR Little Endian transfer syntax (Transfer Syntax UID "1.2.840.10008.1.2"). This means that Pixel Data of Image instances will be retrieved uncompressed. To retrieve instances in any available transfer syntax (typically the one in which instances were originally stored), specify acceptable transfer syntaxes using the wildcard ("application/dicom", "*").

retrieve_instance_frames(study_instance_uid, series_instance_uid, sop_instance_uid, frame_numbers, media_types=None)

Retrieve one or more frames of an image instance.

Parameters

- **study_instance_uid** (str) – Study Instance UID
- **series_instance_uid** (str) – Series Instance UID
- **sop_instance_uid** (str) – SOP Instance UID

- **frame_numbers** (*Sequence[int]*) – One-based positional indices of the frames within the instance
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes

Returns

Pixel data for each frame

Return type

List[bytes]

retrieve_instance_frames_rendered(*study_instance_uid, series_instance_uid, sop_instance_uid, frame_numbers, media_types=None, params=None*)

Retrieve one or more server-side rendered frames of an instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **frame_numbers** (*Sequence[int]*) – One-based positional index of the frame within the instance
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media type (choices: "image/jpeg", "image/jp2", "image/gif", "image/png")
- **params** (*Union[Dict[str, Any], None]*, *optional*) – Additional parameters relevant for given *media_type*, e.g., {"quality": 95} for "image/jpeg" media type

Returns

Rendered representation of frames

Return type

bytes

Note: Not all media types are compatible with all SOP classes.

retrieve_instance_metadata(*study_instance_uid, series_instance_uid, sop_instance_uid*)

Retrieve metadata of an individual instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID

Returns

Metadata of instance in DICOM JSON format

Return type

Dict[str, dict]

retrieve_instance_rendered(*study_instance_uid*, *series_instance_uid*, *sop_instance_uid*,
media_types=None, *params=None*)

Retrieve an individual, server-side rendered instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types (choices: "image/jpeg", "image/jp2", "image/gif", "image/png", "video/gif", "video/mp4", "video/h265", "text/html", "text/plain", "text/xml", "text/rtf", "application/pdf")
- **params** (*Union[Dict[str, Any], None]*, *optional*) – Additional parameters relevant for given *media_type*, e.g., {"quality": 95} for "image/jpeg"

Returns

Rendered representation of instance

Return type

bytes

retrieve_series(*study_instance_uid*, *series_instance_uid*, *media_types=None*)

Retrieve all instances of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the acceptable transfer syntaxes

Returns

Instances

Return type

List[pydicom.dataset.Dataset]

Note: Instances are by default retrieved using Implicit VR Little Endian transfer syntax (Transfer Syntax UID "1.2.840.10008.1.2"). This means that Pixel Data of Image instances will be retrieved uncompressed. To retrieve instances in any available transfer syntax (typically the one in which instances were originally stored), specify acceptable transfer syntaxes using the wildcard ("application/dicom", "*").

retrieve_series_metadata(*study_instance_uid*, *series_instance_uid*)

Retrieve metadata for all instances of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID

Returns

Metadata of instances in DICOM JSON format

Return type

Dict[str, dict]

retrieve_series_rendered(*study_instance_uid*, *series_instance_uid*, *media_types=None*, *params=None*)

Retrieve rendered representation of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types (choices: "image/jpeg", "image/jp2", "image/gif", "image/png", "video/gif", "video/mp4", "video/h265", "text/html", "text/plain", "text/xml", "text/rtf", "application/pdf")
- **params** (*Union[Dict[str, Any], None]*, *optional*) – Additional parameters relevant for given *media_type*, e.g., {"quality": 95} for "image/jpeg"

Returns

Rendered representation of series

Return type

bytes

retrieve_study(*study_instance_uid*, *media_types=None*)

Retrieve all instances of a study.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the acceptable transfer syntaxes

Returns

Instances

Return type

List[pydicom.dataset.Dataset]

Note: Instances are by default retrieved using Implicit VR Little Endian transfer syntax (Transfer Syntax UID "1.2.840.10008.1.2"). This means that Pixel Data of Image instances will be retrieved uncompressed. To retrieve instances in any available transfer syntax (typically the one in which instances were originally stored), specify acceptable transfer syntaxes using the wildcard ("application/dicom", "*").

retrieve_study_metadata(*study_instance_uid*)

Retrieve metadata of all instances of a study.

Parameters**study_instance_uid** (*str*) – Study Instance UID**Returns**

Metadata of instances in DICOM JSON format

Return type

List[Dict[str, dict]]

scheme: `str`

search_for_instances(*study_instance_uid=None, series_instance_uid=None, fuzzymatching=None, limit=None, offset=None, fields=None, search_filters=None, get_remaining=False*)

Search for instances.

Parameters

- **study_instance_uid** (*Union[str, None], optional*) – Study Instance UID
- **series_instance_uid** (*Union[str, None], optional*) – Series Instance UID
- **fuzzymatching** (*Union[bool, None], optional*) – Whether fuzzy semantic matching should be performed
- **limit** (*Union[int, None], optional*) – Maximum number of results that should be returned
- **offset** (*Union[int, None], optional*) – Number of results that should be skipped
- **fields** (*Union[Union[list, tuple, set], None], optional*) – Names of fields (attributes) that should be included in results
- **search_filters** (*Union[Dict[str, Union[str, int, float]], None], optional*) – Search filter criteria as key-value pairs, where *key* is a keyword or a tag of the attribute and *value* is the expected value that should match
- **get_remaining** (*bool, optional*) – Whether remaining results should be included (this may repeatedly query the server for remaining results)

Returns

Instance representations (see [Instance Result Attributes](#))

Return type

List[Dict[str, dict]]

Note: The server may only return a subset of search results. In this case, a warning will notify the client that there are remaining results. Remaining results can be requested via repeated calls using the *offset* parameter.

search_for_series(*study_instance_uid=None, fuzzymatching=None, limit=None, offset=None, fields=None, search_filters=None, get_remaining=False*)

Search for series.

Parameters

- **study_instance_uid** (*Union[str, None], optional*) – Study Instance UID
- **fuzzymatching** (*Union[bool, None], optional*) – Whether fuzzy semantic matching should be performed
- **limit** (*Union[int, None], optional*) – Maximum number of results that should be returned
- **offset** (*Union[int, None], optional*) – Number of results that should be skipped
- **fields** (*Union[Union[list, tuple, set], None], optional*) – Names of fields (attributes) that should be included in results

- **search_filters** (*Union[Dict[str, Union[str, int, float]], None], optional*) – Search filter criteria as key-value pairs, where *key* is a keyword or a tag of the attribute and *value* is the expected value that should match
- **get_remaining** (*bool, optional*) – Whether remaining results should be included (this may repeatedly query the server for remaining results)

Returns

Series representations (see [Series Result Attributes](#))

Return type

List[Dict[str, dict]]

Note: The server may only return a subset of search results. In this case, a warning will notify the client that there are remaining results. Remaining results can be requested via repeated calls using the *offset* parameter.

search_for_studies(*fuzzymatching=None, limit=None, offset=None, fields=None, search_filters=None, get_remaining=False*)

Search for studies.

Parameters

- **fuzzymatching** (*Union[bool, None], optional*) – Whether fuzzy semantic matching should be performed
- **limit** (*Union[int, None], optional*) – Maximum number of results that should be returned
- **offset** (*Union[int, None], optional*) – Number of results that should be skipped
- **fields** (*Union[Sequence[str], None], optional*) – Names of fields (attributes) that should be included in results
- **search_filters** (*Union[dict, None], optional*) – Search filter criteria as key-value pairs, where *key* is a keyword or a tag of the attribute and *value* is the expected value that should match
- **get_remaining** (*bool, optional*) – Whether remaining results should be included (this may repeatedly query the server for remaining results)

Returns

Study representations (see [Study Result Attributes](#))

Return type

List[Dict[str, dict]]

Note: The server may only return a subset of search results. In this case, a warning will notify the client that there are remaining results. Remaining results can be requested via repeated calls using the *offset* parameter.

store_instances(*datasets, study_instance_uid=None*)

Store instances.

Parameters

- **datasets** (*Sequence[pydicom.dataset.Dataset]*) – Instances that should be stored
- **study_instance_uid** (*Union[str, None], optional*) – Study Instance UID

Returns

Information about status of stored instances

Return type

pydicom.dataset.Dataset

stow_url_prefix: Optional[str] = None

url_prefix: str

wado_url_prefix: Optional[str] = None

```
class dicomweb_client.api.DICOMfileClient(url, update_db=False, recreate_db=False, in_memory=False,
                                          db_dir=None, readonly=False)
```

Bases: object

Client for managing DICOM Part10 files in a DICOMweb-like manner.

Facilitates serverless access to data stored locally on a file system as DICOM Part10 files.

Note: The class internally uses an in-memory database, which is persisted on disk to facilitate faster subsequent data access. However, the implementation details of the database and the structure of any database files stored on the file system may change at any time and should not be relied on.

Note: This is **not** an implementation of the DICOM File Service and does not depend on the presence of DICOMDIR files.

base_url

Unique resource locator of the DICOMweb service

Type

str

scheme

Name of the scheme (file)

Type

str

protocol

Name of the protocol (None)

Type

str

url_prefix

URL path prefix for DICOMweb services (part of *base_url*)

Type

str

qido_url_prefix

URL path prefix for QIDO-RS (not part of *base_url*)

Type

Union[str, None]

wado_url_prefix

URL path prefix for WADO-RS (not part of *base_url*)

Type

Union[str, None]

stow_url_prefix

URL path prefix for STOW-RS (not part of *base_url*)

Type

Union[str, None]

delete_url_prefix

URL path prefix for DELETE (not part of *base_url*)

Type

Union[str, None]

Instantiate client.

Parameters

- **url** (*str*) – Unique resource locator of directory where data is stored (must have file scheme)
- **update_db** (*bool, optional*) – Whether the database should be updated (default: False). If True, the client will search *base_dir* recursively for new DICOM Part10 files and create database entries for each file. The client will further delete any database entries for files that no longer exist on the file system.
- **recreate_db** (*bool, optional*) – Whether the database should be recreated (default: False). If True, the client will search *base_dir* recursively for DICOM Part10 files and create database entries for each file.
- **in_memory** (*bool, optional*) – Whether the database should only be stored in memory (default: False).
- **db_dir** (*Union[pathlib.Path, str, None], optional*) – Path to directory where database files should be stored (defaults to *base_dir*)
- **readonly** (*bool, optional*) – Whether data should be considered read-only. Attempts to store or delete data will be denied.

delete_instance(*study_instance_uid, series_instance_uid, sop_instance_uid*)

Delete specified instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID

Return type

None

delete_series(*study_instance_uid, series_instance_uid*)

Delete all instances of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID

- **series_instance_uid** (*str*) – Series Instance UID

Return type

None

delete_study(*study_instance_uid*)

Delete all instances of a study.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID

Return type

None

iter_bulkdata(*url*, *media_types=None*, *byte_range=None*)

Iterate over bulk data items at a given location.

Parameters

- **url** (*str*) – Location of the bulk data
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes
- **byte_range** (*Union[Tuple[int, int], None]*, *optional*) – Start and end of byte range

Returns

Bulk data items

Return type

Iterator[bytes]

Raises**IOError** – When requested resource is not found at *url***iter_instance_frames**(*study_instance_uid*, *series_instance_uid*, *sop_instance_uid*, *frame_numbers*, *media_types=None*)

Iterate over frames of an image instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **frame_numbers** (*List[int]*) – Frame numbers
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes

Returns

Frames

Return type

Iterator[bytes]

iter_series(*study_instance_uid*, *series_instance_uid*, *media_types=None*)

Iterate over all instances of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes

Returns

Instances

Return type

Iterator[pydicom.dataset.Dataset]

iter_study(*study_instance_uid*, *media_types=None*)

Iterate over all instances of a study.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes

Returns

Instances

Return type

Iterator[pydicom.dataset.Dataset]

static lookup_keyword(*tag*)

Look up the keyword of a DICOM attribute.

Parameters**tag** (*Union[str, int, Tuple[int, int], pydicom.tag.BaseTag]*) – Attribute tag (e.g. "00080018")**Returns**

Attribute keyword (e.g. "SOPInstanceUID")

Return type

str

static lookup_tag(*keyword*)

Look up the tag of a DICOM attribute.

Parameters**keyword** (*str*) – Attribute keyword (e.g. "SOPInstanceUID")**Returns**

Attribute tag as HEX string (e.g. "00080018")

Return type

str

retrieve_bulkdata(*url*, *media_types=None*, *byte_range=None*)

Retrieve bulk data at a given location.

Parameters

- **url** (*str*) – Location of the bulk data

- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes
- **byte_range** (*Union[Tuple[int, int], None]*, *optional*) – Start and end of byte range

Returns

Bulk data items

Return type

Iterator[bytes]

Raises

IOError – When requested resource is not found at *url*

retrieve_instance(*study_instance_uid, series_instance_uid, sop_instance_uid, media_types=None*)

Retrieve metadata of a single instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes

Returns

Instance

Return type

pydicom.dataset.Dataset

retrieve_instance_frames(*study_instance_uid, series_instance_uid, sop_instance_uid, frame_numbers, media_types=None*)

Retrieve one or more frames of an image instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **frame_numbers** (*List[int]*) – Frame numbers
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes

Returns

Frames

Return type

List[bytes]

retrieve_instance_frames_rendered(*study_instance_uid, series_instance_uid, sop_instance_uid, frame_numbers, media_types=None, params=None*)

Retrieve server-side rendered frames of an image instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **frame_numbers** (*List[int]*) – Frame numbers
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None], optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes
- **params** (*Union[Dict[str, str], None], optional*) – Additional query parameters

Returns

Rendered representation of frames

Return type

bytes

retrieve_instance_metadata(*study_instance_uid, series_instance_uid, sop_instance_uid*)

Retrieve metadata of a single instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID

Returns

Metadata of instance

Return type

Dict[str, Any]

retrieve_instance_rendered(*study_instance_uid, series_instance_uid, sop_instance_uid, media_types=None, params=None*)

Retrieve an individual, server-side rendered instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None], optional*) – Acceptable media types (choices: "image/jpeg", "image/jp2", "image/gif", "image/png", "video/gif", "video/mp4", "video/h265", "text/html", "text/plain", "text/xml", "text/rtf", "application/pdf")
- **params** (*Union[Dict[str, Any], None], optional*) – Additional parameters relevant for given *media_type*, e.g., {"quality": 95} for "image/jpeg"

Returns

Rendered representation of instance

Return type

bytes

Note: Only rendering of single-frame image instances is currently supported.

retrieve_series(*study_instance_uid*, *series_instance_uid*, *media_types=None*)

Retrieve all instances of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes

Returns

Instances

Return type

Sequence[pydicom.dataset.Dataset]

retrieve_series_metadata(*study_instance_uid*, *series_instance_uid*)

Retrieve metadata of instances in a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID

Returns

Metadata of each instance in series

Return type

List[Dict[str, Any]]

retrieve_series_rendered(*study_instance_uid*, *series_instance_uid*, *media_types=None*, *params=None*)

Retrieve rendered representation of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types (choices: "image/jpeg", "image/jp2", "image/gif", "image/png", "video/gif", "video/mp4", "video/h265", "text/html", "text/plain", "text/xml", "text/rtf", "application/pdf")
- **params** (*Union[Dict[str, Any], None]*, *optional*) – Additional parameters relevant for given *media_type*, e.g., {"quality": 95} for "image/jpeg"

Returns

Rendered representation of series

Return type

bytes

retrieve_study(*study_instance_uid*, *media_types=None*)

Retrieve all instances of a study.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes

Returns

Instances

Return type

Sequence[pydicom.dataset.Dataset]

retrieve_study_metadata(*study_instance_uid*)

Retrieve metadata of instances in a study.

Parameters**study_instance_uid** (*str*) – Study Instance UID**Returns**

Metadata of each instance in study

Return type

List[Dict[str, Any]]

search_for_instances(*study_instance_uid=None, series_instance_uid=None, fuzzymatching=None, limit=None, offset=None, fields=None, search_filters=None, get_remaining=False*)

Search for instances.

Parameters

- **study_instance_uid** (*Union[str, None]*, *optional*) – Study Instance UID
- **series_instance_uid** (*Union[str, None]*, *optional*) – Series Instance UID
- **fuzzymatching** (*Union[bool, None]*, *optional*) – Whether fuzzy semantic matching should be performed
- **limit** (*Union[int, None]*, *optional*) – Maximum number of results that should be returned
- **offset** (*Union[int, None]*, *optional*) – Number of results that should be skipped
- **fields** (*Union[Sequence[str], None]*, *optional*) – Names of fields (attributes) that should be included in results
- **search_filters** (*Union[dict, None]*, *optional*) – Search filter criteria as key-value pairs, where *key* is a keyword or a tag of the attribute and *value* is the expected value that should match
- **get_remaining** (*bool*, *optional*) – Whether remaining results should be included

ReturnsInstances (see [Instance Result Attributes](#))**Return type**

List[Dict[str, dict]]

Note: No additional *fields* are currently supported.

search_for_series(*study_instance_uid=None, fuzzymatching=None, limit=None, offset=None, fields=None, search_filters=None, get_remaining=False*)

Search for series.

Parameters

- **study_instance_uid** (*Union[str, None], optional*) – Study Instance UID
- **fuzzymatching** (*Union[bool, None], optional*) – Whether fuzzy semantic matching should be performed
- **limit** (*Union[int, None], optional*) – Maximum number of results that should be returned
- **offset** (*Union[int, None], optional*) – Number of results that should be skipped
- **fields** (*Union[Sequence[str], None], optional*) – Names of fields (attributes) that should be included in results
- **search_filters** (*Union[dict, None], optional*) – Search filter criteria as key-value pairs, where *key* is a keyword or a tag of the attribute and *value* is the expected value that should match
- **get_remaining** (*bool, optional*) – Whether remaining results should be included

Returns

Series (see [Series Result Attributes](#))

Return type

List[Dict[str, dict]]

search_for_studies(*fuzzymatching=None, limit=None, offset=None, fields=None, search_filters=None, get_remaining=False*)

Search for studies.

Parameters

- **fuzzymatching** (*Union[bool, None], optional*) – Whether fuzzy semantic matching should be performed
- **limit** (*Union[int, None], optional*) – Maximum number of results that should be returned
- **offset** (*Union[int, None], optional*) – Number of results that should be skipped
- **fields** (*Union[Sequence[str], None], optional*) – Names of fields (attributes) that should be included in results
- **search_filters** (*Union[dict, None], optional*) – Search filter criteria as key-value pairs, where *key* is a keyword or a tag of the attribute and *value* is the expected value that should match
- **get_remaining** (*bool, optional*) – Whether remaining results should be included

Returns

Studies (see [Study Result Attributes](#))

Return type

List[Dict[str, dict]]

Note: No additional *fields* are currently supported.

store_instances(*datasets*, *study_instance_uid=None*)

Store instances.

Parameters

- **datasets** (*Sequence[pydicom.dataset.Dataset]*) – Instances that should be stored
- **study_instance_uid** (*Union[str, None]*, *optional*) – Study Instance UID

Returns

Information about status of stored instances

Return type

pydicom.dataset.Dataset

class dicomweb_client.api.DICOMwebClient(*url*, *session=None*, *qido_url_prefix=None*,
wado_url_prefix=None, *stow_url_prefix=None*,
delete_url_prefix=None, *proxies=None*, *headers=None*,
callback=None, *chunk_size=1000000*)

Bases: object

Class for connecting to and interacting with a DICOMweb RESTful service.

base_url

Unique resource locator of the DICOMweb service

Type

str

scheme

Name of the scheme, e.g. "https"

Type

str

protocol

Name of the protocol, e.g. "https"

Type

str

host

IP address or DNS name of the machine that hosts the server

Type

str

port

Number of the port to which the server listens

Type

int

url_prefix

URL path prefix for DICOMweb services (part of *base_url*)

Type

str

qido_url_prefix

URL path prefix for QIDO-RS (not part of *base_url*)

Type

Union[str, None]

wado_url_prefixURL path prefix for WADO-RS (not part of *base_url*)**Type**

Union[str, None]

stow_url_prefixURL path prefix for STOW-RS (not part of *base_url*)**Type**

Union[str, None]

delete_url_prefixURL path prefix for DELETE (not part of *base_url*)**Type**

Union[str, None]

Instantiate client.

Parameters

- **url** (*str*) – Unique resource locator of the DICOMweb service consisting of protocol, host-name (IP address or DNS name) of the machine that hosts the service and optionally port number and path prefix
- **session** (*Union[requests.Session, None], optional*) – Session required to make connections to the DICOMweb service (see `dicomweb_client.session_utils` module to create a valid session if necessary)
- **qido_url_prefix** (*Union[str, None], optional*) – URL path prefix for QIDO RESTful services
- **wado_url_prefix** (*Union[str, None], optional*) – URL path prefix for WADO RESTful services
- **stow_url_prefix** (*Union[str, None], optional*) – URL path prefix for STOW RESTful services
- **delete_url_prefix** (*Union[str, None], optional*) – URL path prefix for DELETE RESTful services
- **proxies** (*Union[Dict[str, str], None], optional*) – Mapping of protocol or protocol + host to the URL of a proxy server
- **headers** (*Union[Dict[str, str], None], optional*) – Custom headers that should be included in request messages, e.g., authentication tokens
- **callback** (*Union[Callable[[requests.Response, ...], requests.Response], None], optional*) – Callback function to manipulate responses generated from requests (see [requests event hooks](#))
- **chunk_size** (*int, optional*) – Maximum number of bytes that should be transferred per data chunk when streaming data from the server using chunked transfer encoding (used by `iter_*`(`)` methods as well as the `store_instances()` method); defaults to $10^{*}6$ bytes (1MB)

Warning: Modifies the passed *session* (in particular header fields), so be careful when reusing the session outside the scope of an instance.

Warning: Choose the value of *chunk_size* carefully. A small value may cause significant network communication and message parsing overhead.

delete_instance(*study_instance_uid*, *series_instance_uid*, *sop_instance_uid*)

Delete specified instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID

Note: The Delete Instance resource is not part of the DICOM standard and may not be supported by all origin servers.

Warning: This method performs a DELETE and should be used with caution.

Return type

None

delete_series(*study_instance_uid*, *series_instance_uid*)

Delete all instances of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID

Note: The Delete Series resource is not part of the DICOM standard and may not be supported by all origin servers.

Warning: This method performs a DELETE and should be used with caution.

Return type

None

delete_study(*study_instance_uid*)

Delete all instances of a study.

Parameters

study_instance_uid (*str*) – Study Instance UID

Note: The Delete Study resource is not part of the DICOM standard and may not be supported by all origin servers.

Warning: This method performs a DELETE and should be used with caution.

Return type

None

iter_bulkdata(*url*, *media_types=None*, *byte_range=None*)

Iterate over bulk data items at a given location.

Parameters

- **url** (*str*) – Location of the bulk data
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes
- **byte_range** (*Union[Tuple[int, int], None]*, *optional*) – Start and end of byte range

Returns

Bulk data items

Return type

Iterator[bytes]

Note: Data is streamed from the DICOMweb server.

iter_instance_frames(*study_instance_uid*, *series_instance_uid*, *sop_instance_uid*, *frame_numbers*, *media_types=None*)

Iterate over frames of an image instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **frame_numbers** (*Sequence[int]*) – One-based positional indices of the frames within the instance
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes

Returns

Pixel data for each frame

Return type

Iterator[bytes]

Note: Data is streamed from the DICOMweb server.

iter_series(*study_instance_uid*, *series_instance_uid*, *media_types=None*)

Iterate over all instances of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the acceptable transfer syntaxes

Returns

Instances

Return type

Iterator[pydicom.dataset.Dataset]

Note: Instances are by default retrieved using Implicit VR Little Endian transfer syntax (Transfer Syntax UID "1.2.840.10008.1.2"). This means that Pixel Data of Image instances will be retrieved uncompressed. To retrieve instances in any available transfer syntax (typically the one in which instances were originally stored), specify acceptable transfer syntaxes using the wildcard ("application/dicom", "*").

Note: Data is streamed from the DICOMweb server.

iter_study(*study_instance_uid*, *media_types=None*)

Iterate over all instances of a study.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the acceptable transfer syntaxes

Returns

Instances

Return type

Iterator[pydicom.dataset.Dataset]

Note: Instances are by default retrieved using Implicit VR Little Endian transfer syntax (Transfer Syntax UID "1.2.840.10008.1.2"). This means that Pixel Data of Image instances will be retrieved uncompressed. To retrieve instances in any available transfer syntax (typically the one in which instances were originally stored), specify acceptable transfer syntaxes using the wildcard ("application/dicom", "*").

Note: Data is streamed from the DICOMweb server.

static lookup_keyword(*tag*)

Look up the keyword of a DICOM attribute.

Parameters

tag (*Union[str, int, Tuple[str, str], pydicom.tag.BaseTag]*) – Attribute tag (e.g. "00080018")

Returns

Attribute keyword (e.g. "SOPInstanceUID")

Return type

str

static lookup_tag(*keyword*)

Look up the tag of a DICOM attribute.

Parameters

keyword (*str*) – Attribute keyword (e.g. "SOPInstanceUID")

Returns

Attribute tag as HEX string (e.g. "00080018")

Return type

str

retrieve_bulkdata(*url, media_types=None, byte_range=None*)

Retrieve bulk data at a given location.

Parameters

- **url** (*str*) – Location of the bulk data
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None], optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes
- **byte_range** (*Union[Tuple[int, int], None], optional*) – Start and end of byte range

Returns

Bulk data items

Return type

Iterator[bytes]

retrieve_instance(*study_instance_uid, series_instance_uid, sop_instance_uid, media_types=None*)

Retrieve an individual instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None], optional*) – Acceptable media types and optionally the UIDs of the acceptable transfer syntaxes

Returns

Instance

Return type

pydicom.dataset.Dataset

Note: Instances are by default retrieved using Implicit VR Little Endian transfer syntax (Transfer Syntax UID "1.2.840.10008.1.2"). This means that Pixel Data of Image instances will be retrieved uncompressed. To retrieve instances in any available transfer syntax (typically the one in which instances were originally stored), specify acceptable transfer syntaxes using the wildcard ("application/dicom", "*").

retrieve_instance_frames(*study_instance_uid*, *series_instance_uid*, *sop_instance_uid*, *frame_numbers*, *media_types=None*)

Retrieve one or more frames of an image instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **frame_numbers** (*Sequence[int]*) – One-based positional indices of the frames within the instance
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the corresponding transfer syntaxes

Returns

Pixel data for each frame

Return type

List[bytes]

retrieve_instance_frames_rendered(*study_instance_uid*, *series_instance_uid*, *sop_instance_uid*, *frame_numbers*, *media_types=None*, *params=None*)

Retrieve one or more server-side rendered frames of an instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **frame_numbers** (*Sequence[int]*) – One-based positional index of the frame within the instance
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media type (choices: "image/jpeg", "image/jp2", "image/gif", "image/png")
- **params** (*Union[Dict[str, Any], None]*, *optional*) – Additional parameters relevant for given *media_type*, e.g., {"quality": 95} for "image/jpeg" media type

Returns

Rendered representation of frames

Return type

bytes

Note: Not all media types are compatible with all SOP classes.

retrieve_instance_metadata(*study_instance_uid*, *series_instance_uid*, *sop_instance_uid*)

Retrieve metadata of an individual instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID

Returns

Metadata of instance in DICOM JSON format

Return type

Dict[str, dict]

retrieve_instance_rendered(*study_instance_uid*, *series_instance_uid*, *sop_instance_uid*,
media_types=None, *params=None*)

Retrieve an individual, server-side rendered instance.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **sop_instance_uid** (*str*) – SOP Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types (choices: "image/jpeg", "image/jp2", "image/gif", "image/png", "video/gif", "video/mp4", "video/h265", "text/html", "text/plain", "text/xml", "text/rtf", "application/pdf")
- **params** (*Union[Dict[str, Any], None]*, *optional*) – Additional parameters relevant for given *media_type*, e.g., {"quality": 95} for "image/jpeg"

Returns

Rendered representation of instance

Return type

bytes

retrieve_series(*study_instance_uid*, *series_instance_uid*, *media_types=None*)

Retrieve all instances of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the acceptable transfer syntaxes

Returns

Instances

Return type

List[pydicom.dataset.Dataset]

Note: Instances are by default retrieved using Implicit VR Little Endian transfer syntax (Transfer Syntax UID "1.2.840.10008.1.2"). This means that Pixel Data of Image instances will be retrieved uncompressed. To retrieve instances in any available transfer syntax (typically the one in which instances were originally stored), specify acceptable transfer syntaxes using the wildcard ("application/dicom", "*").

retrieve_series_metadata(*study_instance_uid*, *series_instance_uid*)

Retrieve metadata for all instances of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID

Returns

Metadata of instances in DICOM JSON format

Return type

Dict[str, dict]

retrieve_series_rendered(*study_instance_uid*, *series_instance_uid*, *media_types=None*, *params=None*)

Retrieve rendered representation of a series.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **series_instance_uid** (*str*) – Series Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types (choices: "image/jpeg", "image/jp2", "image/gif", "image/png", "video/gif", "video/mp4", "video/h265", "text/html", "text/plain", "text/xml", "text/rtf", "application/pdf")
- **params** (*Union[Dict[str, Any], None]*, *optional*) – Additional parameters relevant for given *media_type*, e.g., {"quality": 95} for "image/jpeg"

Returns

Rendered representation of series

Return type

bytes

retrieve_study(*study_instance_uid*, *media_types=None*)

Retrieve all instances of a study.

Parameters

- **study_instance_uid** (*str*) – Study Instance UID
- **media_types** (*Union[Tuple[Union[str, Tuple[str, str]], ...], None]*, *optional*) – Acceptable media types and optionally the UIDs of the acceptable transfer syntaxes

Returns

Instances

Return type

List[pydicom.dataset.Dataset]

Note: Instances are by default retrieved using Implicit VR Little Endian transfer syntax (Transfer Syntax UID "1.2.840.10008.1.2"). This means that Pixel Data of Image instances will be retrieved uncompressed. To retrieve instances in any available transfer syntax (typically the one in which instances were originally stored), specify acceptable transfer syntaxes using the wildcard ("application/dicom", "*").

retrieve_study_metadata(*study_instance_uid*)

Retrieve metadata of all instances of a study.

Parameters**study_instance_uid** (*str*) – Study Instance UID**Returns**

Metadata of instances in DICOM JSON format

Return type

List[Dict[str, dict]]

search_for_instances(*study_instance_uid=None, series_instance_uid=None, fuzzymatching=None, limit=None, offset=None, fields=None, search_filters=None, get_remaining=False*)

Search for instances.

Parameters

- **study_instance_uid** (*Union[str, None], optional*) – Study Instance UID
- **series_instance_uid** (*Union[str, None], optional*) – Series Instance UID
- **fuzzymatching** (*Union[bool, None], optional*) – Whether fuzzy semantic matching should be performed
- **limit** (*Union[int, None], optional*) – Maximum number of results that should be returned
- **offset** (*Union[int, None], optional*) – Number of results that should be skipped
- **fields** (*Union[Union[list, tuple, set], None], optional*) – Names of fields (attributes) that should be included in results
- **search_filters** (*Union[Dict[str, Union[str, int, float]], None], optional*) – Search filter criteria as key-value pairs, where *key* is a keyword or a tag of the attribute and *value* is the expected value that should match
- **get_remaining** (*bool, optional*) – Whether remaining results should be included (this may repeatedly query the server for remaining results)

ReturnsInstance representations (see [Instance Result Attributes](#))**Return type**

List[Dict[str, dict]]

Note: The server may only return a subset of search results. In this case, a warning will notify the client that there are remaining results. Remaining results can be requested via repeated calls using the *offset*

parameter.

search_for_series(*study_instance_uid=None, fuzzymatching=None, limit=None, offset=None, fields=None, search_filters=None, get_remaining=False*)

Search for series.

Parameters

- **study_instance_uid** (*Union[str, None], optional*) – Study Instance UID
- **fuzzymatching** (*Union[bool, None], optional*) – Whether fuzzy semantic matching should be performed
- **limit** (*Union[int, None], optional*) – Maximum number of results that should be returned
- **offset** (*Union[int, None], optional*) – Number of results that should be skipped
- **fields** (*Union[Union[list, tuple, set], None], optional*) – Names of fields (attributes) that should be included in results
- **search_filters** (*Union[Dict[str, Union[str, int, float]], None], optional*) – Search filter criteria as key-value pairs, where *key* is a keyword or a tag of the attribute and *value* is the expected value that should match
- **get_remaining** (*bool, optional*) – Whether remaining results should be included (this may repeatedly query the server for remaining results)

Returns

Series representations (see [Series Result Attributes](#))

Return type

List[Dict[str, dict]]

Note: The server may only return a subset of search results. In this case, a warning will notify the client that there are remaining results. Remaining results can be requested via repeated calls using the *offset* parameter.

search_for_studies(*fuzzymatching=None, limit=None, offset=None, fields=None, search_filters=None, get_remaining=False*)

Search for studies.

Parameters

- **fuzzymatching** (*Union[bool, None], optional*) – Whether fuzzy semantic matching should be performed
- **limit** (*Union[int, None], optional*) – Maximum number of results that should be returned
- **offset** (*Union[int, None], optional*) – Number of results that should be skipped
- **fields** (*Union[Sequence[str], None], optional*) – Names of fields (attributes) that should be included in results
- **search_filters** (*Union[dict, None], optional*) – Search filter criteria as key-value pairs, where *key* is a keyword or a tag of the attribute and *value* is the expected value that should match

- **get_remaining** (*bool, optional*) – Whether remaining results should be included (this may repeatedly query the server for remaining results)

Returns

Study representations (see [Study Result Attributes](#))

Return type

List[Dict[str, dict]]

Note: The server may only return a subset of search results. In this case, a warning will notify the client that there are remaining results. Remaining results can be requested via repeated calls using the *offset* parameter.

set_chunk_size(*chunk_size*)

Set value of *chunk_size* attribute.

Parameters

chunk_size (*int*) – Maximum number of bytes that should be transferred per data chunk when streaming data from the server using chunked transfer encoding (used by *iter_**() methods as well as the *store_instances()* method)

Return type

None

set_http_retry_params(*retry=True, max_attempts=5, wait_exponential_multiplier=1000, retrievable_error_codes=(<HTTPStatus.TOO_MANY_REQUESTS: 429>, <HTTPStatus.REQUEST_TIMEOUT: 408>, <HTTPStatus.SERVICE_UNAVAILABLE: 503>, <HTTPStatus.GATEWAY_TIMEOUT: 504>)*)

Set parameters for HTTP retrying logic.

These parameters determine whether and how individual HTTP requests will be retried in case the origin server responds with an error code defined in *retrievable_error_codes*. The retrying method uses exponential back off using the multiplier *wait_exponential_multiplier* for a max attempts defined by *max_attempts*.

Parameters

- **retry** (*bool, optional*) – Whether HTTP retrying should be performed, if it is set to *False*, the rest of the parameters are ignored.
- **max_attempts** (*int, optional*) – The maximum number of request attempts.
- **wait_exponential_multiplier** (*float, optional*) – Exponential multiplier applied to delay between attempts in ms.
- **retrievable_error_codes** (*tuple, optional*) – Tuple of HTTP error codes to retry if raised.

Return type

None

store_instances(*datasets, study_instance_uid=None*)

Store instances.

Parameters

- **datasets** (*Sequence[pydicom.dataset.Dataset]*) – Instances that should be stored
- **study_instance_uid** (*Union[str, None], optional*) – Study Instance UID

Returns

Information about status of stored instances

Return type

pydicom.dataset.Dataset

7.1.2 dicomweb_client.cli module

Command Line Interface (CLI)

`dicomweb_client.cli.main(args)`

Main entry point for the `dicomweb_client` command line program.

`dicomweb_client`

Client for DICOMweb RESTful services.

```
usage: dicomweb_client [-h] [-v] [-u NAME] [-p PASSWORD] [--ca CERT-FILE]
                    [--cert CERT-FILE] [--bearer-token TOKEN] [--url URL]
                    [--chunk-size NUM]
                    {search,retrieve,store} ...
```

-h, --help

show this help message and exit

-v, --verbosity

logging verbosity that maps to a logging level (default: error, -v: warning, -vv: info, -vvv: debug, -vvvv: debug + traceback); all log messages are written to standard error

-u <name>, --user <name>

username for authentication with the DICOMweb service

-p <password>, --password <password>

password for authentication with the DICOMweb service

--ca <cert-file>

path to a CA bundle file

--cert <cert-file>

path to a client certificate file in PEM format

--bearer-token <token>

bearer token for authentication with the DICOMweb service

--url <url>

uniform resource locator of the DICOMweb service

--chunk-size <num>

maximum size of a network transfer chunk in bytes

dicomweb_client retrieve

WADO-RS: Web Access to DICOM Objects by RESTful Services.

```
usage: dicomweb_client retrieve [-h] INFORMATION ENTITIES ...
```

-h, --help

show this help message and exit

dicomweb_client retrieve bulkdata

Retrieve bulk data of a DICOM object from a known location.

```
usage: dicomweb_client retrieve bulkdata [-h]
                                     [--media-type MEDIATYPE [MEDIATYPE ...]]
                                     --uri URI
```

-h, --help

show this help message and exit

--media-type <mediatype>

acceptable media type and the optionally the UID of a corresponding transfer syntax separated by a whitespace (e.g., “image/jpeg” or “image/jpeg 1.2.840.10008.1.2.4.50”)

--uri <uri>

unique resource identifier of bulk data element

dicomweb_client retrieve instances

Retrieve data for an individual DICOM instance.

```
usage: dicomweb_client retrieve instances [-h] --study UID --series UID
                                     --instance UID
                                     {metadata,full,frames} ...
```

-h, --help

show this help message and exit

--study <uid>

unique study identifier (StudyInstanceUID)

--series <uid>

unique series identifier (SeriesInstanceUID)

--instance <uid>

unique instance identifier (SOPInstanceUID)

dicomweb_client retrieve instances frames

Retrieve one or more frames of the pixel data element of an individual DICOM instance.

```
usage: dicomweb_client retrieve instances frames [-h] [--save]
                                             [--output-dir PATH]
                                             [--media-type MEDIATYPE [MEDIATYPE ...]]
                                             [--numbers NUM [NUM ...]]
                                             [--show]
```

-h, --help

show this help message and exit

--save

whether downloaded data should be saved

--output-dir <path>

path to directory where downloaded data should be saved

--media-type <mediatype>

acceptable media type and the optionally the UID of a corresponding transfer syntax separated by a whitespace (e.g., “image/jpeg” or “image/jpeg 1.2.840.10008.1.2.4.50”)

--numbers <num>

frame numbers

--show

display retrieved images

dicomweb_client retrieve instances full

Retrieve a DICOM instance.

```
usage: dicomweb_client retrieve instances full [-h] [--save]
                                             [--output-dir PATH]
                                             [--media-type MEDIATYPE [MEDIATYPE ...]]
```

-h, --help

show this help message and exit

--save

whether downloaded data should be saved

--output-dir <path>

path to directory where downloaded data should be saved

--media-type <mediatype>

acceptable media type and the optionally the UID of a corresponding transfer syntax separated by a whitespace (e.g., “image/jpeg” or “image/jpeg 1.2.840.10008.1.2.4.50”)

dicomweb_client retrieve instances metadata

Retrieve metadata of an individual DICOM instance.

```
usage: dicomweb_client retrieve instances metadata [-h]
                                                [--prettyfy | --dicomize]
                                                [--save]
                                                [--output-dir PATH]
```

-h, --help

show this help message and exit

--prettyfy

pretty print JSON response message

--dicomize

convert JSON response message to DICOM data set

--save

whether downloaded data should be saved

--output-dir <path>

path to directory where downloaded data should be saved

dicomweb_client retrieve series

Retrieve data for all DICOM instances of a given DICOM series.

```
usage: dicomweb_client retrieve series [-h] --study UID --series UID
                                       {metadata,full} ...
```

-h, --help

show this help message and exit

--study <uid>

unique study identifier (StudyInstanceUID)

--series <uid>

unique series identifier (SeriesInstanceUID)

dicomweb_client retrieve series full

Retrieve DICOM instances of a given DICOM series.

```
usage: dicomweb_client retrieve series full [-h] [--save] [--output-dir PATH]
                                             [--media-type MEDIATYPE [MEDIATYPE ...]]
```

-h, --help

show this help message and exit

--save

whether downloaded data should be saved

--output-dir <path>

path to directory where downloaded data should be saved

--media-type <mediatype>

acceptable media type and the optionally the UID of a corresponding transfer syntax separated by a whitespace (e.g., “image/jpeg” or “image/jpeg 1.2.840.10008.1.2.4.50”)

dicomweb_client retrieve series metadata

Retrieve metadata of DICOM instances of a given DICOM series.

```
usage: dicomweb_client retrieve series metadata [-h] [--prettyfy | --dicomize]
        [--save] [--output-dir PATH]
```

-h, --help

show this help message and exit

--prettyfy

pretty print JSON response message

--dicomize

convert JSON response message to DICOM data set

--save

whether downloaded data should be saved

--output-dir <path>

path to directory where downloaded data should be saved

dicomweb_client retrieve studies

Retrieve data for all DICOM instances of a given DICOM study.

```
usage: dicomweb_client retrieve studies [-h] --study UID {metadata,full} ...
```

-h, --help

show this help message and exit

--study <uid>

unique study identifier (StudyInstanceUID)

dicomweb_client retrieve studies full

Retrieve DICOM instances of a given DICOM study.

```
usage: dicomweb_client retrieve studies full [-h] [--save] [--output-dir PATH]
        [--media-type MEDIATYPE [MEDIATYPE ...]]
```

-h, --help

show this help message and exit

--save

whether downloaded data should be saved

--output-dir <path>

path to directory where downloaded data should be saved

--media-type <mediatype>

acceptable media type and the optionally the UID of a corresponding transfer syntax separated by a whitespace (e.g., "image/jpeg" or "image/jpeg 1.2.840.10008.1.2.4.50")

dicomweb_client retrieve studies metadata

Retrieve metadata of DICOM instances of a given DICOM study.

```
usage: dicomweb_client retrieve studies metadata [-h]
                                                [--prettyify | --dicomize]
                                                [--save] [--output-dir PATH]
```

-h, --help

show this help message and exit

--prettyify

pretty print JSON response message

--dicomize

convert JSON response message to DICOM data set

--save

whether downloaded data should be saved

--output-dir <path>

path to directory where downloaded data should be saved

dicomweb_client search

QIDO-RS: Query based on ID for DICOM Objects by RESTful Services.

```
usage: dicomweb_client search [-h] INFORMATION ENTITIES ...
```

-h, --help

show this help message and exit

dicomweb_client search instances

Search for DICOM instances.

```
usage: dicomweb_client search instances [-h] [--prettyify | --dicomize]
                                         [--filter KEY=VALUE] [--field NAME]
                                         [--limit NUM] [--offset NUM] [--fuzzy]
                                         [--study UID] [--series UID]
```

- h, --help**
show this help message and exit
- pretty**
pretty print JSON response message
- dicomize**
convert JSON response message to DICOM data set
- filter <key=value>**
query filter criterion
- field <name>**
field that should be included in response
- limit <num>**
number of items that should be maximally retrieved
- offset <num>**
number of items that should be skipped
- fuzzy**
perform fuzzy matching
- study <uid>**
unique study identifier (StudyInstanceUID)
- series <uid>**
unique series identifier (SeriesInstanceUID)

dicomweb_client search series

Search for DICOM series.

```
usage: dicomweb_client search series [-h] [--filter KEY=VALUE] [--field NAME]
                                     [--limit NUM] [--offset NUM] [--fuzzy]
                                     [--pretty | --dicomize] [--study UID]
```

- h, --help**
show this help message and exit
- filter <key=value>**
query filter criterion
- field <name>**
field that should be included in response
- limit <num>**
number of items that should be maximally retrieved
- offset <num>**
number of items that should be skipped
- fuzzy**
perform fuzzy matching

--prettyfy

pretty print JSON response message

--dicomize

convert JSON response message to DICOM data set

--study <uid>

unique study identifier (StudyInstanceUID)

dicomweb_client search studies

Search for DICOM studies.

```
usage: dicomweb_client search studies [-h] [--filter KEY=VALUE] [--field NAME]
                                     [--limit NUM] [--offset NUM] [--fuzzy]
                                     [--prettyfy | --dicomize]
```

-h, --help

show this help message and exit

--filter <key=value>

query filter criterion

--field <name>

field that should be included in response

--limit <num>

number of items that should be maximally retrieved

--offset <num>

number of items that should be skipped

--fuzzy

perform fuzzy matching

--prettyfy

pretty print JSON response message

--dicomize

convert JSON response message to DICOM data set

dicomweb_client store

STOW-RS: Store Over the Web by RESTful Services.

```
usage: dicomweb_client store [-h] INFORMATION ENTITIES ...
```

-h, --help

show this help message and exit

dicomweb_client store instances

Store DICOM instances.

```
usage: dicomweb_client store instances [-h] [--study UID] PATH [PATH ...]
```

path

paths to DICOM files that should be loaded

-h, --help

show this help message and exit

--study <uid>

unique study identifier (StudyInstanceUID)

7.1.3 dicomweb_client.log module

Utility functions for logging configuration

`dicomweb_client.log.configure_logging(verbosity)`

Configures the root logger with a “stderr” stream handler that directs logging messages to standard error (to allow capturing program standard output, e.g. in order to redirect it to a file).

Logging verbosity maps to levels as follows:

```
0 -> no messages
1 -> CRITICAL, ERROR & WARN/WARNING messages
2 -> CRITICAL, ERROR, WARN/WARNING, & INFO messages
3 -> CRITICAL, ERROR, WARN/WARNING, INFO & DEBUG messages
4 -> all messages
```

Parameters

verbosity (*int*) – logging verbosity

Returns

package root logger

Return type

`logging.Logger`

7.1.4 dicomweb_client.session_utils module

`dicomweb_client.session_utils.add_certs_to_session(session, ca_bundle=None, cert=None)`

Adds CA bundle and certificate to an existing session.

Parameters

- **session** (*requests.Session*) – input session
- **ca_bundle** (*str, optional*) – path to CA bundle file
- **cert** (*str, optional*) – path to client certificate file in Privacy Enhanced Mail (PEM) format

Returns

verified session

Return type

requests.Session

`dicomweb_client.session_utils.create_session()`

Creates an unauthorized session.

Returns

unauthorized session

Return type

requests.Session

`dicomweb_client.session_utils.create_session_from_auth(auth)`

Creates a session from a given AuthBase object.

Parameters

auth (*requests.auth.AuthBase*) – an implementation of *requests.auth.AuthBase* to be used for authentication with services

Returns

authorized session

Return type

requests.Session

`dicomweb_client.session_utils.create_session_from_gcp_credentials(google_credentials=None)`

Creates an authorized session for Google Cloud Platform.

Parameters

google_credentials (*Any*) – Google cloud credentials. (see <https://cloud.google.com/docs/authentication/production> for more information on Google cloud authentication). If not set, will be initialized to `google.auth.default()`

Returns

Google cloud authorized session

Return type

requests.Session

`dicomweb_client.session_utils.create_session_from_user_pass(username, password)`

Creates a session from a given username and password.

Parameters

- **username** (*str*) – username for authentication with services
- **password** (*str*) – password for authentication with services

Returns

authorized session

Return type

requests.Session

7.1.5 dicomweb_client.uri module

Utilities for DICOMweb URI manipulation.

```
class dicomweb_client.uri.URI(base_url, study_instance_uid=None, series_instance_uid=None,  
sop_instance_uid=None, frames=None, suffix=None, permissive=False)
```

Bases: object

Class to represent a fully qualified HTTP[S] URI to a DICOMweb resource.

<http://dicom.nema.org/medical/dicom/current/output/html/part18.html>

This is an immutable class. Use *URI.update()* to create copies of an instance with updated (new) values for its attributes.

Given an HTTP[S] *base_url*, a valid DICOMweb-compatible URI would be:

- <base_url> (no DICOMweb suffix)
- <base_url>/studies/<study_instance_uid>
- <base_url>/studies/<study_instance_uid>/metadata
- <base_url>/studies/<study_instance_uid>/rendered
- <base_url>/studies/<study_instance_uid>/thumbnail
- <base_url>/studies/<study_instance_uid>/series/<series_instance_uid>
- <base_url>/studies/<study_instance_uid>/series/<series_instance_uid>/metadata
- <base_url>/studies/<study_instance_uid>/series/<series_instance_uid>/rendered
- <base_url>/studies/<study_instance_uid>/series/<series_instance_uid>/thumbnail
- <base_url>/studies/<study_instance_uid>/series/<series_instance_uid>/instances/
<sop_instance_uid>
- <base_url>/studies/<study_instance_uid>/series/<series_instance_uid>/instances/
<sop_instance_uid>/metadata
- <base_url>/studies/<study_instance_uid>/series/<series_instance_uid>/instances/
<sop_instance_uid>/rendered
- <base_url>/studies/<study_instance_uid>/series/<series_instance_uid>/instances/
<sop_instance_uid>/thumbnail
- <base_url>/studies/<study_instance_uid>/series/<series_instance_uid>/instances/
<sop_instance_uid>/frames/<frames>
- <base_url>/studies/<study_instance_uid>/series/<series_instance_uid>/instances/
<sop_instance_uid>/frames/<frames>/rendered
- <base_url>/studies/<study_instance_uid>/series/<series_instance_uid>/instances/
<sop_instance_uid>/frames/<frames>/thumbnail

Instantiates an object.

As per the DICOM Standard, the Study, Series, and Instance UIDs must be a series of numeric components (0-9) separated by the period . character, with a maximum length of 64 characters. If the *permissive* flag is set to True, any alpha-numeric or special characters (except for / and @) may be used.

Parameters

- **base_url** (*str*) – DICOMweb service HTTP[S] URL. Trailing forward slashes are not permitted.

- **study_instance_uid** (*str, optional*) – DICOM Study Instance UID.
- **series_instance_uid** (*str, optional*) – DICOM Series Instance UID.
- **sop_instance_uid** (*str, optional*) – DICOM SOP Instance UID.
- **frames** (*Sequence[int], optional*) – A non-empty sequence of positive frame numbers in ascending order.
- **suffix** (*URISuffix, optional*) – Suffix attached to the DICOM resource URI. This could refer to a metadata, rendered, or thumbnail resource.
- **permissive** (*bool*) – If True, relaxes the DICOM Standard validation for UIDs (see main docstring for details). This option is made available since users may be occasionally forced to work with DICOMs or services that may be in violation of the standard. Unless required, use of this flag is **not** recommended, since non-conformant UIDs may lead to unexpected errors downstream, e.g., rejection by a DICOMweb server, etc.

Raises

ValueError – In the following cases:

- *base_url* has a trailing slash. - *base_url* does not use the HTTP[S] addressing scheme. - *base_url* is incompatible with the DICOMweb standard. - *series_instance_uid* is supplied without *study_instance_uid*. - *sop_instance_uid* is supplied without *study_instance_uid* or *series_instance_uid*. - *frames* is supplied without *study_instance_uid*, *series_instance_uid*, or *sop_instance_uid*. - *frames* is empty. - A frame number in *frames* is not positive. - The frame numbers in *frames* are not in ascending order. - *suffix* is *URISuffix.METADATA* with *frames* or without *study_instance_uid*. - *suffix* is *URISuffix.RENDERED* without *study_instance_uid*. - *suffix* is *URISuffix.THUMBNAIL* without *study_instance_uid*. - Any one of *study_instance_uid*, *series_instance_uid*, or *sop_instance_uid* does not meet the DICOM Standard UID spec in the docstring.

base_uri()

Returns *URI* for the DICOM Service within this object.

Return type

dicomweb_client.uri.URI

property base_url: str

Returns the Base (DICOMweb Service) URL.

Return type

str

frame_uri()

Returns *URI* for the DICOM frames within this object.

Return type

dicomweb_client.uri.URI

property frames: Optional[Tuple[int, ...]]

Returns the sequence of frame numbers, if available.

Return type

typing.Optional[typing.Tuple[int, ...]]

classmethod from_string(dicomweb_uri, uri_type=None, permissive=False)

Parses the string to return the URI.

Any valid DICOMweb compatible HTTP[S] URI is permitted, e.g., <SERVICE>/studies/<StudyInstanceUID>/series/<SeriesInstanceUID>.

Parameters

- **dicomweb_uri** (*str*) – An HTTP[S] DICOMweb-compatible URI.
- **uri_type** (*URIType*, *optional*) – The expected DICOM resource type referenced by the object. If set, it validates that the resource-scope of the *dicomweb_uri* matches the expected type.
- **permissive** (*bool*) – Set if permissive handling of UIDs (if any) in *dicomweb_uri* is required. See the class initializer docstring for details.

Returns

The newly constructed *URI* object.

Return type

URI

Raises

ValueError – If the URI cannot be parsed or the actual URI type doesn't match the specified expected *uri_type*.

instance_uri()

Returns *URI* for the DICOM Instances within this object.

Return type

dicomweb_client.uri.URI

property parent: *URI*

Returns a URI to the “parent” resource.

Depending on the *type* of the current *URI*, the *URI* of the parent resource is defined as:

Current	Suffixed	Parent
Service	N/A	Service
Study	No	Service
Study	Yes	Study
Series	No	Study
Series	Yes	Series
Instance	No	Series
Instance	Yes	Instance
Frame	No	Instance
Frame	Yes	Frame

Returns

An instance pointing to the parent resource.

Return type

URI

property permissive: *bool*

Returns the permissive parameter value in the initializer.

Return type

bool

property series_instance_uid: *Optional[str]*

Returns the Series UID, if available.

Return type

typing.Optional[str]

series_uri()Returns *URI* for the DICOM Series within this object.**Return type***dicomweb_client.uri.URI***property sop_instance_uid: Optional[str]**

Returns the Instance UID, if available.

Return type

typing.Optional[str]

property study_instance_uid: Optional[str]

Returns the Study UID, if available.

Return type

typing.Optional[str]

study_uri()Returns *URI* for the DICOM Study within this object.**Return type***dicomweb_client.uri.URI***property suffix: Optional[URISuffix]**

Returns the DICOM resource suffix, if available.

Return typetyping.Optional[*dicomweb_client.uri.URISuffix*]**property type: URIType**The *URIType* of DICOM resource referenced by the object.**Return type***dicomweb_client.uri.URIType***update**(*base_url=None, study_instance_uid=None, series_instance_uid=None, sop_instance_uid=None, frames=None, suffix=None, permissive=False*)Creates a new *URI* object based on the current one.Replaces the specified *URI* components in the current *URI* to create the new one.**Parameters**

- **base_url** (*str, optional*) – DICOMweb service HTTP[S] URL to use in the new *URI* or *None* if the *base_url* from the current *URI* should be used.
- **study_instance_uid** (*str, optional*) – Study Instance UID to use in the new *URI* or *None* if the *study_instance_uid* from the current *URI* should be used.
- **series_instance_uid** (*str, optional*) – Series Instance UID to use in the new *URI* or *None* if the *series_instance_uid* from the current *URI* should be used.
- **sop_instance_uid** (*str, optional*) – SOP Instance UID to use in the new *URI* or *None* if the *sop_instance_uid* from the current *URI* should be used.
- **frames** (*Sequence[int], optional*) – Frame numbers to use in the new *URI* or *None* if the *frames* from the current *URI* should be used.

- **suffix** (*URISuffix*, *optional*) – Suffix to use in the new *URI* or *None* if the *suffix* from the current *URI* should be used.
- **permissive** (*bool*, *optional*) – Set if permissive handling of UIDs (if any) in the updated *URI* is required. See the class initializer docstring for details.

Returns

The newly constructed *URI* object.

Return type

URI

Raises

ValueError – If the new *URI* is invalid (e.g., if only the SOP Instance UID is specified, but the Series Instance UID is missing in the current *URI*).

```
class dicomweb_client.uri.URISuffix(value)
```

Bases: Enum

Optional suffixes for a DICOM resource.

```
METADATA = 'metadata'
```

```
RENDERED = 'rendered'
```

```
THUMBNAIL = 'thumbnail'
```

```
class dicomweb_client.uri.URIType(value)
```

Bases: Enum

Type of DICOM resource the *URI* points to.

```
FRAME = 'frame'
```

```
INSTANCE = 'instance'
```

```
SERIES = 'series'
```

```
SERVICE = 'service'
```

```
STUDY = 'study'
```

```
dicomweb_client.uri.build_query_string(params=None)
```

Build query string for a request message.

Parameters

params (*Union[Dict[str, Any], None]*, *optional*) – Query parameters as mapping of key-value pairs; in case a key should be included more than once with different values, values need to be provided in form of an iterable (e.g., {"key": ["value1", "value2"]} will result in "?key=value1&key=value2")

Returns

Query string

Return type

str

```
dicomweb_client.uri.parse_query_parameters(fuzzymatching=None, limit=None, offset=None,
                                           fields=None, search_filters=None)
```

Parse query parameters for inclusion into a query string.

Parameters

- **fuzzymatching** (*Union[bool, None], optional*) – Whether fuzzy semantic matching should be performed
- **limit** (*Union[int, None], optional*) – Maximum number of results that should be returned
- **offset** (*Union[int, None], optional*) – Number of results that should be skipped
- **fields** (*Union[Sequence[str], None], optional*) – Names of fields (attributes) that should be included in results
- **search_filters** (*Union[Dict[str, Any], None], optional*) – Search filter criteria as key-value pairs, where *key* is a keyword or a tag of the attribute and *value* is the expected value that should match

Returns

Sanitized and sorted query parameters

Return type

collections.OrderedDict

7.1.6 dicomweb_client.ext.gcp.session_utils module

7.1.7 dicomweb_client.ext.gcp.uri module

Utilities for Google Cloud Healthcare DICOMweb API URI manipulation.

For details, visit: <https://cloud.google.com/healthcare>

class `dicomweb_client.ext.gcp.uri.GoogleCloudHealthcareURL`(*project_id, location, dataset_id, dicom_store_id*)

Bases: object

URL container for DICOM Stores under the [Google Cloud Healthcare API](#).

This class facilitates the parsing and creation of `URI.base_url` corresponding to DICOMweb API Service URLs under the `v1` API. The URLs are of the form: `https://healthcare.googleapis.com/v1/projects/{project_id}/locations/{location}/datasets/{dataset_id}/dicomStores/{dicom_store_id}/dicomWeb`

project_id

The ID of the [GCP Project](#) that contains the DICOM Store.

Type

str

location

The [Region name](#) of the geographic location configured for the Dataset that contains the DICOM Store.

Type

str

dataset_id

The ID of the [Dataset](#) that contains the DICOM Store.

Type

str

dicom_store_id

The ID of the [DICOM Store](#).

Type

str

dataset_id: str

dicom_store_id: str

classmethod from_string(*base_url*)

Create an instance from *base_url*.

Parameters

base_url (*str*) – The URL for the DICOMweb API Service endpoint corresponding to a [CHC API DICOM Store](#). See class docstring for supported formats.

Raises

ValueError – If *base_url* does not match the specifications in the class docstring.

Return type

dicomweb_client.ext.gcp.uri.GoogleCloudHealthcareURL

location: str

project_id: str

INDICES AND TABLES

- genindex
- modindex
- search

PYTHON MODULE INDEX

d

`dicomweb_client.api`, 23
`dicomweb_client.cli`, 53
`dicomweb_client.ext.gcp.uri`, 68
`dicomweb_client.log`, 61
`dicomweb_client.session_utils`, 61
`dicomweb_client.uri`, 63

Symbols

- bearer-token
 - dicomweb_client command line option, 53
- ca
 - dicomweb_client command line option, 53
- cert
 - dicomweb_client command line option, 53
- chunk-size
 - dicomweb_client command line option, 53
- dicomize
 - dicomweb_client-retrieve-instances-metadata command line option, 56
 - dicomweb_client-retrieve-series-metadata command line option, 57
 - dicomweb_client-retrieve-studies-metadata command line option, 58
 - dicomweb_client-search-instances command line option, 59
 - dicomweb_client-search-series command line option, 60
 - dicomweb_client-search-studies command line option, 60
- field
 - dicomweb_client-search-instances command line option, 59
 - dicomweb_client-search-series command line option, 59
 - dicomweb_client-search-studies command line option, 60
- filter
 - dicomweb_client-search-instances command line option, 59
 - dicomweb_client-search-series command line option, 59
 - dicomweb_client-search-studies command line option, 60
- fuzzy
 - dicomweb_client-search-instances command line option, 59
 - dicomweb_client-search-series command line option, 59
 - dicomweb_client-search-studies command line option, 60
- line option, 60
- help
 - dicomweb_client command line option, 53
 - dicomweb_client-retrieve command line option, 54
 - dicomweb_client-retrieve-bulkdata command line option, 54
 - dicomweb_client-retrieve-instances command line option, 54
 - dicomweb_client-retrieve-instances-frames command line option, 55
 - dicomweb_client-retrieve-instances-full command line option, 55
 - dicomweb_client-retrieve-instances-metadata command line option, 56
 - dicomweb_client-retrieve-series command line option, 56
 - dicomweb_client-retrieve-series-full command line option, 56
 - dicomweb_client-retrieve-series-metadata command line option, 57
 - dicomweb_client-retrieve-studies command line option, 57
 - dicomweb_client-retrieve-studies-full command line option, 57
 - dicomweb_client-retrieve-studies-metadata command line option, 58
 - dicomweb_client-search command line option, 58
 - dicomweb_client-search-instances command line option, 58
 - dicomweb_client-search-series command line option, 59
 - dicomweb_client-search-studies command line option, 60
 - dicomweb_client-store command line option, 60
 - dicomweb_client-store-instances command line option, 61
- instance
 - dicomweb_client-retrieve-instances command line option, 54

```

--limit
  dicomweb_client-search-instances
    command line option, 59
  dicomweb_client-search-series command
    line option, 59
  dicomweb_client-search-studies command
    line option, 60
--media-type
  dicomweb_client-retrieve-bulkdata
    command line option, 54
  dicomweb_client-retrieve-instances-frames
    command line option, 55
  dicomweb_client-retrieve-instances-full
    command line option, 55
  dicomweb_client-retrieve-series-full
    command line option, 57
  dicomweb_client-retrieve-studies-full
    command line option, 58
--numbers
  dicomweb_client-retrieve-instances-frames
    command line option, 55
--offset
  dicomweb_client-search-instances
    command line option, 59
  dicomweb_client-search-series command
    line option, 59
  dicomweb_client-search-studies command
    line option, 60
--output-dir
  dicomweb_client-retrieve-instances-frames
    command line option, 55
  dicomweb_client-retrieve-instances-full
    command line option, 55
  dicomweb_client-retrieve-instances-metadata
    command line option, 56
  dicomweb_client-retrieve-series-full
    command line option, 56
  dicomweb_client-retrieve-series-metadata
    command line option, 57
  dicomweb_client-retrieve-studies-full
    command line option, 58
  dicomweb_client-retrieve-studies-metadata
    command line option, 58
--password
  dicomweb_client command line option, 53
--prettify
  dicomweb_client-retrieve-instances-metadata
    command line option, 56
  dicomweb_client-retrieve-series-metadata
    command line option, 57
  dicomweb_client-retrieve-studies-metadata
    command line option, 58
  dicomweb_client-search-instances
    command line option, 59
  dicomweb_client-search-series command
    line option, 59
  dicomweb_client-search-studies command
    line option, 60
--save
  dicomweb_client-retrieve-instances-frames
    command line option, 55
  dicomweb_client-retrieve-instances-full
    command line option, 55
  dicomweb_client-retrieve-instances-metadata
    command line option, 56
  dicomweb_client-retrieve-series-full
    command line option, 56
  dicomweb_client-retrieve-series-metadata
    command line option, 57
  dicomweb_client-retrieve-studies-full
    command line option, 57
  dicomweb_client-retrieve-studies-metadata
    command line option, 58
--series
  dicomweb_client-retrieve-instances
    command line option, 54
  dicomweb_client-retrieve-series command
    line option, 56
  dicomweb_client-search-instances
    command line option, 59
--show
  dicomweb_client-retrieve-instances-frames
    command line option, 55
--study
  dicomweb_client-retrieve-instances
    command line option, 54
  dicomweb_client-retrieve-series command
    line option, 56
  dicomweb_client-retrieve-studies
    command line option, 57
  dicomweb_client-search-instances
    command line option, 59
  dicomweb_client-search-series command
    line option, 60
  dicomweb_client-store-instances command
    line option, 61
--uri
  dicomweb_client-retrieve-bulkdata
    command line option, 54
--url
  dicomweb_client command line option, 53
--user
  dicomweb_client command line option, 53
--verbosity
  dicomweb_client command line option, 53
-h
  dicomweb_client command line option, 53
  dicomweb_client-retrieve command line

```


option, 54

`dicomweb_client-retrieve-bulkdata` command line option, 54

`dicomweb_client-retrieve-instances` command line option, 54

`dicomweb_client-retrieve-instances-frames` command line option, 55

`dicomweb_client-retrieve-instances-full` command line option, 55

`dicomweb_client-retrieve-instances-metadata` command line option, 56

`dicomweb_client-retrieve-series` command line option, 56

`dicomweb_client-retrieve-series-full` command line option, 56

`dicomweb_client-retrieve-series-metadata` command line option, 57

`dicomweb_client-retrieve-studies` command line option, 57

`dicomweb_client-retrieve-studies-full` command line option, 57

`dicomweb_client-retrieve-studies-metadata` command line option, 58

`dicomweb_client-search` command line option, 58

`dicomweb_client-search-instances` command line option, 58

`dicomweb_client-search-series` command line option, 59

`dicomweb_client-search-studies` command line option, 60

`dicomweb_client-store` command line option, 60

`dicomweb_client-store-instances` command line option, 61

`-p dicomweb_client` command line option, 53

`-u dicomweb_client` command line option, 53

`-v dicomweb_client` command line option, 53

A

`add_certs_to_session()` (in module `dicomweb_client.session_utils`), 61

B

`base_uri()` (`dicomweb_client.uri.URI` method), 64

`base_url` (`dicomweb_client.api.DICOMClient` attribute), 23

`base_url` (`dicomweb_client.api.DICOMfileClient` attribute), 32

`base_url` (`dicomweb_client.api.DICOMwebClient` attribute), 41

`base_url` (`dicomweb_client.uri.URI` property), 64

`build_query_string()` (in module `dicomweb_client.uri`), 67

C

`configure_logging()` (in module `dicomweb_client.log`), 61

`create_session()` (in module `dicomweb_client.session_utils`), 62

`create_session_from_auth()` (in module `dicomweb_client.session_utils`), 62

`create_session_from_gcp_credentials()` (in module `dicomweb_client.session_utils`), 62

`create_session_from_user_pass()` (in module `dicomweb_client.session_utils`), 62

D

`dataset_id` (`dicomweb_client.ext.gcp.uri.GoogleCloudHealthcareURL` attribute), 68, 69

`delete_instance()` (`dicomweb_client.api.DICOMClient` method), 23

`delete_instance()` (`dicomweb_client.api.DICOMfileClient` method), 33

`delete_instance()` (`dicomweb_client.api.DICOMwebClient` method), 43

`delete_series()` (`dicomweb_client.api.DICOMClient` method), 23

`delete_series()` (`dicomweb_client.api.DICOMfileClient` method), 33

`delete_series()` (`dicomweb_client.api.DICOMwebClient` method), 43

`delete_study()` (`dicomweb_client.api.DICOMClient` method), 24

`delete_study()` (`dicomweb_client.api.DICOMfileClient` method), 34

`delete_study()` (`dicomweb_client.api.DICOMwebClient` method), 43

`delete_url_prefix` (`dicomweb_client.api.DICOMClient` attribute), 24

`delete_url_prefix` (`dicomweb_client.api.DICOMfileClient` attribute), 33

`delete_url_prefix` (`dicomweb_client.api.DICOMwebClient` attribute), 42

`dicom_store_id` (`dicomweb_client.ext.gcp.uri.GoogleCloudHealthcareURL` attribute), 68, 69

`DICOMClient` (class in `dicomweb_client.api`), 23

DICOMfileClient (*class in dicomweb_client.api*), 32

dicomweb_client command line option

- bearer-token, 53
- ca, 53
- cert, 53
- chunk-size, 53
- help, 53
- password, 53
- url, 53
- user, 53
- verbosity, 53
- h, 53
- p, 53
- u, 53
- v, 53

dicomweb_client.api module, 23

dicomweb_client.cli module, 53

dicomweb_client.ext.gcp.uri module, 68

dicomweb_client.log module, 61

dicomweb_client.session_utils module, 61

dicomweb_client.uri module, 63

dicomweb_client-retrieve command line option

- help, 54
- h, 54

dicomweb_client-retrieve-bulkdata command line option

- help, 54
- media-type, 54
- uri, 54
- h, 54

dicomweb_client-retrieve-instances command line option

- help, 54
- instance, 54
- series, 54
- study, 54
- h, 54

dicomweb_client-retrieve-instances-frames command line option

- help, 55
- media-type, 55
- numbers, 55
- output-dir, 55
- save, 55
- show, 55
- h, 55

dicomweb_client-retrieve-instances-full command line option

- help, 55
- media-type, 55
- output-dir, 55
- save, 55
- h, 55

dicomweb_client-retrieve-instances-metadata command line option

- dicomize, 56
- help, 56
- output-dir, 56
- prettify, 56
- save, 56
- h, 56

dicomweb_client-retrieve-series command line option

- help, 56
- series, 56
- study, 56
- h, 56

dicomweb_client-retrieve-series-full command line option

- help, 56
- media-type, 57
- output-dir, 56
- save, 56
- h, 56

dicomweb_client-retrieve-series-metadata command line option

- dicomize, 57
- help, 57
- output-dir, 57
- prettify, 57
- save, 57
- h, 57

dicomweb_client-retrieve-studies command line option

- help, 57
- study, 57
- h, 57

dicomweb_client-retrieve-studies-full command line option

- help, 57
- media-type, 58
- output-dir, 58
- save, 57
- h, 57

dicomweb_client-retrieve-studies-metadata command line option

- dicomize, 58
- help, 58
- output-dir, 58
- prettify, 58

--save, 58
 -h, 58
 dicomweb_client-search command line option
 --help, 58
 -h, 58
 dicomweb_client-search-instances command
 line option
 --dicomize, 59
 --field, 59
 --filter, 59
 --fuzzy, 59
 --help, 58
 --limit, 59
 --offset, 59
 --prettify, 59
 --series, 59
 --study, 59
 -h, 58
 dicomweb_client-search-series command line
 option
 --dicomize, 60
 --field, 59
 --filter, 59
 --fuzzy, 59
 --help, 59
 --limit, 59
 --offset, 59
 --prettify, 59
 --study, 60
 -h, 59
 dicomweb_client-search-studies command line
 option
 --dicomize, 60
 --field, 60
 --filter, 60
 --fuzzy, 60
 --help, 60
 --limit, 60
 --offset, 60
 --prettify, 60
 -h, 60
 dicomweb_client-store command line option
 --help, 60
 -h, 60
 dicomweb_client-store-instances command
 line option
 --help, 61
 --study, 61
 -h, 61
 path, 61
 DICOMwebClient (class in dicomweb_client.api), 41

F

FRAME (dicomweb_client.uri.URIType attribute), 67

frame_uri() (dicomweb_client.uri.URI method), 64
 frames (dicomweb_client.uri.URI property), 64
 from_string() (dicomweb_client.ext.gcp.uri.GoogleCloudHealthcareURL
 class method), 69
 from_string() (dicomweb_client.uri.URI class
 method), 64

G

GoogleCloudHealthcareURL (class in di-
 comweb_client.ext.gcp.uri), 68

H

host (dicomweb_client.api.DICOMwebClient attribute),
 41

I

INSTANCE (dicomweb_client.uri.URIType attribute), 67
 instance_uri() (dicomweb_client.uri.URI method), 65
 iter_bulkdata() (dicomweb_client.api.DICOMClient
 method), 24
 iter_bulkdata() (di-
 comweb_client.api.DICOMfileClient method),
 34
 iter_bulkdata() (di-
 comweb_client.api.DICOMwebClient method),
 44
 iter_instance_frames() (di-
 comweb_client.api.DICOMClient method),
 24
 iter_instance_frames() (di-
 comweb_client.api.DICOMfileClient method),
 34
 iter_instance_frames() (di-
 comweb_client.api.DICOMwebClient method),
 44
 iter_series() (dicomweb_client.api.DICOMClient
 method), 25
 iter_series() (dicomweb_client.api.DICOMfileClient
 method), 34
 iter_series() (dicomweb_client.api.DICOMwebClient
 method), 45
 iter_study() (dicomweb_client.api.DICOMClient
 method), 25
 iter_study() (dicomweb_client.api.DICOMfileClient
 method), 35
 iter_study() (dicomweb_client.api.DICOMwebClient
 method), 45

L

location (dicomweb_client.ext.gcp.uri.GoogleCloudHealthcareURL
 attribute), 68, 69

lookup_keyword() (di-
 comweb_client.api.DICOMfileClient static
 method), 35

lookup_keyword() (*dicomweb_client.api.DICOMwebClient static method*), 46

lookup_tag() (*dicomweb_client.api.DICOMfileClient static method*), 35

lookup_tag() (*dicomweb_client.api.DICOMwebClient static method*), 46

M

main() (*in module dicomweb_client.cli*), 53

METADATA (*dicomweb_client.uri.URISuffix attribute*), 67

module

dicomweb_client.api, 23

dicomweb_client.cli, 53

dicomweb_client.ext.gcp.uri, 68

dicomweb_client.log, 61

dicomweb_client.session_utils, 61

dicomweb_client.uri, 63

P

parent (*dicomweb_client.uri.URI property*), 65

parse_query_parameters() (*in module dicomweb_client.uri*), 67

path

dicomweb_client-store-instances command line option, 61

permissive (*dicomweb_client.uri.URI property*), 65

port (*dicomweb_client.api.DICOMwebClient attribute*), 41

project_id (*dicomweb_client.ext.gcp.uri.GoogleCloudHealthcareURL attribute*), 68, 69

protocol (*dicomweb_client.api.DICOMClient attribute*), 25

protocol (*dicomweb_client.api.DICOMfileClient attribute*), 32

protocol (*dicomweb_client.api.DICOMwebClient attribute*), 41

Q

qido_url_prefix (*dicomweb_client.api.DICOMClient attribute*), 26

qido_url_prefix (*dicomweb_client.api.DICOMfileClient attribute*), 32

qido_url_prefix (*dicomweb_client.api.DICOMwebClient attribute*), 41

R

RENDERED (*dicomweb_client.uri.URISuffix attribute*), 67

retrieve_bulkdata() (*dicomweb_client.api.DICOMClient method*), 26

retrieve_bulkdata() (*dicomweb_client.api.DICOMfileClient method*), 35

retrieve_bulkdata() (*dicomweb_client.api.DICOMwebClient method*), 46

retrieve_instance() (*dicomweb_client.api.DICOMClient method*), 26

retrieve_instance() (*dicomweb_client.api.DICOMfileClient method*), 36

retrieve_instance() (*dicomweb_client.api.DICOMwebClient method*), 46

retrieve_instance_frames() (*dicomweb_client.api.DICOMClient method*), 26

retrieve_instance_frames() (*dicomweb_client.api.DICOMfileClient method*), 36

retrieve_instance_frames() (*dicomweb_client.api.DICOMwebClient method*), 47

retrieve_instance_frames_rendered() (*dicomweb_client.api.DICOMClient method*), 27

retrieve_instance_frames_rendered() (*dicomweb_client.api.DICOMfileClient method*), 36

retrieve_instance_frames_rendered() (*dicomweb_client.api.DICOMwebClient method*), 47

retrieve_instance_metadata() (*dicomweb_client.api.DICOMClient method*), 27

retrieve_instance_metadata() (*dicomweb_client.api.DICOMfileClient method*), 37

retrieve_instance_metadata() (*dicomweb_client.api.DICOMwebClient method*), 48

retrieve_instance_rendered() (*dicomweb_client.api.DICOMClient method*), 27

retrieve_instance_rendered() (*dicomweb_client.api.DICOMfileClient method*), 37

retrieve_instance_rendered() (*dicomweb_client.api.DICOMwebClient method*), 48

retrieve_series() (*dicomweb_client.api.DICOMClient method*), 28

retrieve_series()	(dicomweb_client.api.DICOMfileClient method), 38	comweb_client.api.DICOMfileClient method), 39
retrieve_series()	(dicomweb_client.api.DICOMwebClient method), 48	search_for_instances() (dicomweb_client.api.DICOMwebClient method), 50
retrieve_series_metadata()	(dicomweb_client.api.DICOMClient method), 28	search_for_series() (dicomweb_client.api.DICOMClient method), 30
retrieve_series_metadata()	(dicomweb_client.api.DICOMfileClient method), 38	search_for_series() (dicomweb_client.api.DICOMfileClient method), 39
retrieve_series_metadata()	(dicomweb_client.api.DICOMwebClient method), 49	search_for_series() (dicomweb_client.api.DICOMwebClient method), 51
retrieve_series_rendered()	(dicomweb_client.api.DICOMClient method), 29	search_for_studies() (dicomweb_client.api.DICOMClient method), 31
retrieve_series_rendered()	(dicomweb_client.api.DICOMfileClient method), 38	search_for_studies() (dicomweb_client.api.DICOMfileClient method), 40
retrieve_series_rendered()	(dicomweb_client.api.DICOMwebClient method), 49	search_for_studies() (dicomweb_client.api.DICOMwebClient method), 51
retrieve_study()	(dicomweb_client.api.DICOMClient method), 29	SERIES (dicomweb_client.uri.URIType attribute), 67
retrieve_study()	(dicomweb_client.api.DICOMfileClient method), 38	series_instance_uid (dicomweb_client.uri.URI property), 65
retrieve_study()	(dicomweb_client.api.DICOMwebClient method), 49	series_uri() (dicomweb_client.uri.URI method), 66
retrieve_study_metadata()	(dicomweb_client.api.DICOMClient method), 29	SERVICE (dicomweb_client.uri.URIType attribute), 67
retrieve_study_metadata()	(dicomweb_client.api.DICOMfileClient method), 39	set_chunk_size() (dicomweb_client.api.DICOMwebClient method), 52
retrieve_study_metadata()	(dicomweb_client.api.DICOMwebClient method), 50	set_http_retry_params() (dicomweb_client.api.DICOMwebClient method), 52
		sop_instance_uid (dicomweb_client.uri.URI property), 66
		store_instances() (dicomweb_client.api.DICOMClient method), 31
		store_instances() (dicomweb_client.api.DICOMfileClient method), 40
		store_instances() (dicomweb_client.api.DICOMwebClient method), 52
		stow_url_prefix (dicomweb_client.api.DICOMClient attribute), 32
		stow_url_prefix (dicomweb_client.api.DICOMfileClient attribute), 33
		stow_url_prefix (dicomweb_client.api.DICOMwebClient attribute), 42
		STUDY (dicomweb_client.uri.URIType attribute), 67

S

scheme (dicomweb_client.api.DICOMClient attribute), 29
 scheme (dicomweb_client.api.DICOMfileClient attribute), 32
 scheme (dicomweb_client.api.DICOMwebClient attribute), 41
 search_for_instances() (dicomweb_client.api.DICOMClient method), 30
 search_for_instances() (di-

`study_instance_uid` (*dicomweb_client.uri.URI property*), 66
`study_uri()` (*dicomweb_client.uri.URI method*), 66
`suffix` (*dicomweb_client.uri.URI property*), 66

T

`THUMBNAIL` (*dicomweb_client.uri.URISuffix attribute*), 67
`type` (*dicomweb_client.uri.URI property*), 66

U

`update()` (*dicomweb_client.uri.URI method*), 66
`URI` (*class in dicomweb_client.uri*), 63
`URISuffix` (*class in dicomweb_client.uri*), 67
`URIType` (*class in dicomweb_client.uri*), 67
`url_prefix` (*dicomweb_client.api.DICOMClient attribute*), 32
`url_prefix` (*dicomweb_client.api.DICOMfileClient attribute*), 32
`url_prefix` (*dicomweb_client.api.DICOMwebClient attribute*), 41

W

`wado_url_prefix` (*dicomweb_client.api.DICOMClient attribute*), 32
`wado_url_prefix` (*dicomweb_client.api.DICOMfileClient attribute*), 32
`wado_url_prefix` (*dicomweb_client.api.DICOMwebClient attribute*), 42